# SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

SD BIO LINE H.pylori (04FK10)

## English

### About the test

**[Introduction]** *Helicobacter pylori* (*H. pylori*) was initially isolated by Warren and Marshall from biopsy specimens taken from patients suffering from active chronic gastritis. In fact, it is now clear that *H. pylori* is the principle etiologic agent in the pathology of type B gastritis (chronic active antral gastritis), for which it appears to be the triggering and perhaps aggravating factor. Increasing data are being accumulated regarding the fundamental role of *H. pylori* in active chronic gastritis, gastric ulcer and duodenal ulcer and its close correlation with gastric lesions. *H. pylori* is typically isolated in culture medium and examined by microscopy after staining or is detected by a urease test. Both of these techniques take time to implement and their sensitivity and specificity have yet to be demonstrated. Immunochromatographic (rapid) techniques for the detection of antibodies specific to *H. pylori* have substantially resolved these problems, ensuring a serological monitoring in a very short time using simple, highly specific technology without recourse to invasive techniques. The test for *H. pylori* can be utilized as a rapid screening process for large populations of patients and is indicated in early diagnosis of *H. pylori* infection as the immune response can often precede clinical manifestations of disease. From a diagnostic point of view, a high serum level specific antibodies against *H. pylori* should be interpreted as an indication of type B asymptomatic gastritis. **[Principle]** The SD BIOLINE H. pylori test contains a membrane strip, which is pre-coated with *H.pylori* capture antigen on test line region. The *H. pylori* antigen colloid gold conjugate and plasma or serum specimen moves along the membrane chromatographically to the test line (T) and forms a visible line as the antigen-antibody-antigen gold particle complex forms with high degree of sensitivity and specificity. **[Intended Use]** The SD BIOLINE H. pylori test is a rapid test for the qualitative detection of antibodies of all isotypes (IgG, IgM, IgA, etc) specific to *Helicobacter pylori* in human serum or plasma. This test kit is intended as an aid in the diagnosis of H. pylori infection in patients with gastrointestinal symptoms. The SD BIOLINE H. pylori kit is intended only for an initial screening test and reactive specimens should be confirmed by another supplemental assay.

### Materials provided and active ingredients of main components

- SD BIOLINE H. pylori test kit contains following items to perform the assay:
  - 30 Test devices with desiccant in individual foil pouches
  - Assay diluent (1 x 4 mL/vial)
  - 1 Instructions for use
- Active ingredients of main components
  - 1 test device included : Gold Conjugates : *Helicobacter pylori* antigen gold colloid (1.0±0.2 µg). Test line : *Helicobacter pylori* antigen (0.430±0.086 µg), Control line : Goat anti- *Helicobacter pylori* polyclonal antibody (0.65±0.13 µg)
  - Assay diluents includes : 50 mM Tris-HCl buffer (q.s.), Sodium azide (q.s.)

### Materials required but not provided

- Micropipette, Protective gloves, Timer, Biohazard container

### Kit storage and stability

- The SD BIOLINE H.pylori test kit must be stored at a temperature between 2 °C and 30 °C. Do not freeze the kit or its components.
- The test device is sensitive to both heat and humidity.
- Check the humidity indicator on the desiccant for color change and throw the pouch if the color indicates saturation (yellow OK. Discard if green.).
- Perform the test immediately after removing the test device from the foil pouch.
- Do not use the test kit beyond its expiration date.
- The shelf life of the kit is as indicated on the outer package.
- Do not use the test kit if the pouch is damaged or the seal is broken.

### Warnings

- For *in vitro* diagnostic use only. Do not reuse the test device.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens and wash hands thoroughly afterward.
- Avoid splashing or aerosol formation of specimen and assay diluent.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.
- Do not mix or interchange different specimens.
- Do not mix reagent of different lots or those for other products.
- The assay diluent contains a low concentration of sodium azide as a preservative. Sodium azide is toxic and should be handled carefully to avoid ingestion and skin contact.

### Specimen collection and storage

Plasma or Serum

- [Plasma]** Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) and then centrifuge the tube to get plasma specimen.
- [Serum]** Using venipuncture, draw whole blood into the collection tube (NOT containing anticoagulant). Leave blood to coagulate for 30 minutes and then centrifuge the tube to get serum specimen of supernatant.
- If plasma or serum specimens are not tested immediately, they should be refrigerated at 2 - 8 °C. For storage periods longer than 1 week, freezing (below -20 °C) is required. Bring plasma or serum specimen to room temperature (15 - 30 °C) prior to use.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Use standard laboratory techniques to clarify such specimens prior to testing.

Precaution

- Anticoagulants including heparin, EDTA, and sodium citrate do not affect the test result. Use of other anticoagulants have not been validated. Their use may affect the test result.
- Hemolytic specimens, specimens containing rheumatoid factors and lipemic, icteric specimens can interfere with the test and impair results.
- Use a new pipette tip for each specimen in order to avoid cross-contamination of either specimens which could cause erroneous results.

### Test procedure (Refer to figure)

- Bring all kit components and specimens to room temperature between 15 °C and 30 °C prior to opening and testing.
- Remove the test device from the foil pouch and place it on a flat, dry surface.
- Using a micropipette, dispense 10 µl of serum or plasma to the specimen well of the test device, then dispense 3 drops of assay diluent.
- Interpret test results 10 minutes after adding assay diluent.

**Caution :** Do not read the test results after 20 minutes; reading after 20 minutes can yield erroneous results.

### Test interpretation (Refer to figure)

- A purple control line will appear in the left section of the result window to show that the test is working properly.
- The right section of the result window indicates the test results. If another purple line appears in the right section of the result window, this is the test line.
  - Negative result :** The presence of only Control Line ("C") in the result window indicates a negative result.
  - Positive result :** The presence of two colored lines ("T" and "C") in the result window, regardless of which line appears first, indicates a positive result.
- \*Caution:** The presence of any line, no matter how faint, the result is considered positive.
- Invalid result :** No Control Line ("C") in the result window means the test is invalid. The instructions may not have been followed correctly or the test may have been compromised. When an invalid result is obtained, it is recommended that the specimen be retested using a new device.

### Test limitations

- A negative result does not preclude the possibility of infection with *H.pylori*. Other clinically available tests are required if questionable results are obtained.
- As with all diagnostic test, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- The test procedure, precautions and interpretation of results for this test must be followed when testing.

### Internal quality control

The SD BIOLINE H.pylori test device has Test Line ("T") and Control Line ("C") on the surface of the device. Neither the Test Line nor the Control Line is visible in the result window prior to applying a specimen. The Control Line is used for procedural control and show only that the diluent has been applied successfully and that the active ingredients of main components on the strip are functional, but is not a guarantee that the specimen has been properly applied ; it is not a positive specimen control.

### Performance characteristics

- Expected Values**

The majority of individuals exposed to *H. pylori* is developed IgG antibodies against the organism. It is reported that *H. pylori* is universally distributed and an estimated 50% of the world's population is colonized by *H. pylori*. A relatively large proportion of patients who have positive levels of IgG antibodies are asymptomatic even though they are colonized with the *H. pylori*. Therefore, IgG antibodies levels do not necessarily correlate with the severity of clinical symptoms.

- Sensitivity and Specificity**

The clinical trial studies were carried out to evaluate the SD BIOLINE H. pylori.

A total of 122 random serum samples collected from Korea were tested on the SD BIOLINE H.pylori in Korea. In a comparison of the SD BIOLINE H. pylori versus a commercial H. pylori ELISA kit, the BIOLINE H.pylori showed a sensitivity of 95.9% (71/74) and a specificity of 89.6 % (43/48).

Reference assay	SD BIOLINE H.pylori		Total Results
	Positive	Negative	
Commercial H. pylori ELISA	Positive 71	Negative 3	74
Sensitivity (95% CI)	5	43	48
Specificity (95% CI)	95.9 <span> </span> % (88.8 - 98.6 <span> </span> %)	89.6 <span> </span> % (77.8 - 95.5 <span> </span> %)	

- Reproducibility of SD BIOLINE H. pylori test results has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.
- Analytical sensitivity : The limit of detection; the smallest amount of the target marker that can be precisely detected ; have been equal or superior to a leading commercial H.pylori antibody detection rapid test.

<p><i>Product Disclaimer:</i></p> <p><i>While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.</i></p> <p><i>Warning:</i></p> <p><i>The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect positive or negative diagnosis using this product.</i></p>
--

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