# SD BIO LINE Chlamydia

## Explanation of the test

In the western world, *Chlamydia trachomatis* is a bacterium which causes a sexually transmitted infection. Chlamydia is very common disease, which, because of its complications, should be taken very seriously. In women, the commonest site for infection is the cervix, or neck of the womb. It can also affect the urethra, which is the tube which carries urine from the bladder to the outside world. Other sites for chlamydial infection in both men and women include the tissues around the rectum and also the eyes. The most worrying effect of a chlamydial infection in women is that of potential fertility problems, due to inflammation of the fallopian tubes or cervix. The disease is particularly common among young people. One of the most common ways of testing for Chlamydia is for the GP (general practitioner) to collect a cell sample from the infected area (cervix or penis) with a cotton swab, which is sent to a laboratory for evaluation and results by the cell culture methods or direct fluorescent antibody assay, etc. Recently, the immunochromatography format (rapid) to detect Chlamydia antigen are available.

The SD BIOLINE Chlamydia is a solid phase immunochromatographic assay for the rapid, qualitative detection of Chlamydia antigen directly from endocervical swab and cytology brush specimens. This test is intended for professional use as an aid in the presumptive diagnosis of chlamydial infection. To perform the test, an endocervical clinical specimen is obtained and placed into a tube containing Reagent A (Extraction solution), after 2 minutes, Reagents B (Neutralization solution) is added to the tube. After extraction and neutralization, 3 drops of extracted sample is added to the test cassette sample well (S). The SD BIOLINE Chlamydia contains a membrane strip, which is pre-coated with mouse monoclonal anti-*Chlamydia trachomatis* on test band region. The complex of sample including chlamydia antigen and mouse monoclonal anti-*Chlamydia trachomatis*-colloid gold conjugate moves along the membrane chromatographically to the test region and forms a visible line as the antibody-antigen-antibody gold particle complex forms. Therefore, the formation of a visible line in the test region indicates a positive result for the detection of Chlamydia antigen. When the Chlamydia antigen are absent in the sample, no visible color band in the test region.

The SD BIOLINE Chlamydia test device has a letter of T and C as "Test Line" and "Control Line" on the surface of the device. Both the Test Line and Control Line in result window are not visible before applying any samples. The Control Line is used for procedural control. Control line should always appear if test procedure is performed properly and the test reagents of control line are working.

## Materials provided

The SD BIOLINE Chlamydia kit contains following items to perform the assay.

- 25 Test devices individually foil pouched with a desiccant
- Reagent A (1 x 10ml/vial) : Extraction solution
- Reagent B (1 x 20ml/vial) : Neutralization solution
- 50 Sterile swabs, 25 Sample collection tubes
- 25 300µl disposable droppers
- 25 600µl disposable droppers
- 1 Instructions for use

## Precautions / kit storage and stability

- For *in vitro* diagnostic use only. Do not re-use the test device.
- The SD BIOLINE Chlamydia test device should be stored at 2~30 °C. Do not freeze the kit or components.
- Do not use it beyond the expiration.
- To obtain accurate results, you must follow the instruction for use.
- Use appropriate precautions in the collection, handling, storage and disposal of the specimens and used kit contents. Discard used materials in a proper biohazard container.
- Reagent A contains sodium hydroxide - a basic solution, Reagent B contains hydrochloric acid - an acidic solution. If either of the reagents contact the skin or eye, flush with large volumes of water.
- Use only sterile swabs or cytology brushes to obtain endocervical specimens.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Check the desiccant for color change of the humidity indicator and to throw away the pouch if the color indicates saturation.
- The test device is sensitive to humidity and as well as to heat. Perform the test immediately after removing the test device from the pouch.

## Specimen Collection and Storage

- [Swab]** A first swab should be taken to remove excess mucus from the exocervix. The second swab should be inserted in to the endocervical canal, past the squamocolumnar junction, until most of the tip is no longer visible. This will permit acquisition of columnar or cuboidal epithelial cells which are the main reservoir of chlamydia organisms. Firmly rotate the swab for 15 ~ 20 seconds. The swab should be withdrawn without contamination with exocervical or vaginal cells.
- [Cytology brush]** After cleaning the exocervix with the Dacron swab, insert the cytology brush into the endocervical canal past the squamocolumnar junction. Leave in place two to three seconds. Rotate the cytology brush two full turns, withdraw the brush without touching any vaginal surface.  
**Caution :** Do not use cytology brushes with pregnant patients.
- The specimen may be tested immediately or returned to the provided transport tube for storage or transport. Do not place the specimen in any transport device containing medium interferes with the assay and viability of the organisms is not required for the assay.
- Specimens may be stored for 6 hours at 15 ~ 27°C or 72 hours refrigerated 2 ~ 8°C. It is recommended that specimens be processed as soon as possible after collection.

## Procedure of the test

### [Extraction Procedure : preparation of extracted sample]

- Bring all kit components and specimen to room temperature prior to testing.
- Open the empty sample collection tube.
- Take Reagent A up to the Fill Line as shown in Figure (about 300µl). And then, transfer into the tube.
- Insert the patient swab into the tube containing Reagent A (300µl/tube).  
→ Compress the bottom of the tube between the thumb and forefinger and twirl the swab 10 times.
- Wait 2 minutes.  
→ Compress the bottom of the tube between the thumb and forefinger and twirl the swab 10 times.
- Hold the dropper vertically, draw reagent B solution up to the Fill Line as shown in Figure. (about 600µl) And then, transfer 600µl of reagent B into the tube.  
→ Compress the bottom of the tube between the thumb and forefinger and twirl the swab 10 times.
- Express the liquid from the swab by compressing the middle of the tube and pulling the swab up through it. Discard the swab.

### [Test Procedure]

- Bring all kit components and specimen to room temperature prior to testing.
- Remove the test device from the foil pouch and place it on a clean, dry surface.
- Assemble dropping cap on the sample collection tube, add 3 drops (about 110µl) of the extracted sample from the tube to the round Sample well (S) on the test device.

- As the test begins to work, you will see purple color move across the membrane in the center of the test device.
- Interpret test results at 15 minutes. Some positive results may be seen earlier.

### Caution

- The above interpreting time is based on reading the test results at room temperature (15~30°C). If your room temperature is significantly lower than 15°C, then the interpreting time should be properly increased.
- All reagents, test device and clinical specimens should be at room temperature before beginning the assay.

## Interpretation of the test

**Negative Result :** The presence of only one purple color band ("Control Line") within the result window indicates a negative result. A negative result indicates that the specimen is presumptive negative for the presence of chlamydial antigen.

**Positive Result :** The presence of two color bands ("Test Line" and "Control Line") within the result window, no matter which band appears first, indicates a positive result. A positive result indicates that the specimen is presumptive positive for the presence of chlamydial antigen.

**Invalid Result :** If the purple color band ("Control Line") is not visible within the result window after performing the test, the result is considered invalid. The directions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be re-tested.

## Limitations of the test

- The SD BIOLINE Chlamydia test has been tested using endocervical swab and cytology brush clinical specimens for the qualitative detection of chlamydia antigen. Performance with other specimens has not been assessed.
- Detection of chlamydia is dependent on the number of organisms present in the specimen. This may be affected by specimen collection methods and patient factors such as age, history of STD, presence of symptoms, etc. The minimum detection level of this test may vary according to serovar.
- Test results should be interpreted in conjunction with other laboratory and clinical data available to the physician. Standard chlamydial cell culture methods should be used in the evaluation of suspected sexual abuse and the other medicolegal cases where diagnosis could lead to adverse psychosocial impact.

## Internal Quality Control

The SD BIOLINE Chlamydia test device has a letter of T and C as "Test Line" and "Control Line" on the surface of the case. Both the Test Line and Control Line in result window are not visible before applying any samples. The Control Line is used for procedural control. The control line of the RDT only shows that the diluents has been applied successfully, and that the active ingredients of main components on the strip was still functional, but is not a guarantee that the sample has been properly applied and does not represent a positive sample control.

## Performance Characteristics

- Comparison with culture method

The accuracy of SD BIOLINE Chlamydia test was evaluated in comparison to culture positive/negative results of endocervical specimens. Of the 71 culture positive specimens, SD BIOLINE Chlamydia test correctly identified 93.1% (67/72) and of the 647 culture negative specimens, SD BIOLINE Chlamydia test correctly identified 98.8% (639/647)

Reference	SD BIOLINE Chlamydia			Total Results
	Results	Positive	Negative	
culture	Positive	67	5	72
	Negative	8	639	647
<b>Total Results</b>		75	644	719

- Comparison study with a commercially available Rapid Chlamydia test

The accuracy of SD BIOLINE Chlamydia test was also evaluated against a commercially available rapid Chlamydia test using endocervical specimens. The used specimens were tested and confirmed with culture and DFA method as positive (51)/negative(123). Of the 51 positive specimens and 123 negative specimens, SD BIOLINE Chlamydia test correctly identified 100%.

Reference	SD BIOLINE Chlamydia			Total Results
	Results	Positive	Negative	
USA product	Positive	51	0	51
	Negative	0	123	123
<b>Total Results</b>		51	123	174

## Reproducibility Studies

Reproducibility of SD BIOLINE Chlamydia has been demonstrated by studies (within-run, between-run and batch-to-batch) with in-house reference panels. All values were identical to reference panel acceptance criteria.

## Bibliography of suggested reading

- Marjorie J. Miller, The Laboratory Diagnosis of Chlamydia trachomatis and Genital Mycoplasmas. Journal of Medical Technology, 2:8 August, 1985. p508-512
- International Application Published under the Patent Cooperation Treaty (PCT), International Publication Number : WO89/08262, 8 September 1989, World Intellectual Property Organization.

Date Issued : 2013. 10  
09FK10-01-En-0

### Disclaimer:

Whilst every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the Manufacturer and Distributor and the result may accordingly be affected by environmental factors and / or user error. A person who is the subject of the diagnosis should consult a doctor for further confirmation of the result.

### Warning:

The Manufacturers and Distributors of this product shall not be liable for any losses, liability, claims, costs or damages whether direct or indirect or consequential arising out of or related to an incorrect diagnosis, whether positive or negative, in the use of this product

Manufactured by

**STANDARD DIAGNOSTICS, INC.**


65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do, Republic of Korea  
Tel : 82-31-899-2800 Fax : 82-31-899-2840  
http://www.standardia.com sales@standardia.com

# One Step Chlamydia Antigen Rapid Test


# SD BIO LINE Chlamydia

## Test Procedure


**1** Add 300µl of Reagent A


**2** Extraction


**3** Add 600µl of Reagent B


**4** Neutralization


**5** Discard the swab.

**6** Assemble dropping cap on the sample collection tube.


**7** Add 3 drops (about 110µl) of extracted sample.

## Interpretation


15 min

## How to interpret test results


### Negative

One line "C" in the result window.



### Positive

Two lines "C", "T" in the result window.

### Invalid

No "C" line in the result window.  
It is recommended that the specimen should be re-tested.

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