REF 62FK10

Panbio™
COVID-19/Flu A&B
Rapid Panel
(Nasopharyngeal)

Abbott Rapid Diagnostics Jena GmbH
Orlaweg 1, D-07743 Jena, Germany
www.globalpointofcare.abbott

Date Issued: 2022.04
62FK10A-08-A0

© 2022 Abbott. All rights reserved.

All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Scan the QR Code for more information about this product. / Escanee el código QR para obtener más información sobre este producto. / Escaneie o código QR para mais informações deste produto. / Scanați codul QR pentru mai multe informații despre acest produs. / Skenirajte QR kod za više informacija o ovom proizvodu. / Для получения дополнительной информации об этом продукте отсканируйте QR-код. / Zeskanuj kod QR, aby uzyskać więcej informacji o produkcie. / Další informace o tomto výrobku získáte naskenováním QR kódu.

Introduction

The coronavirus disease (COVID-19) is an infectious disease caused by a newly discovered coronavirus, severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2).¹ The SARS-CoV-2 is a β -coronavirus, which is an enveloped non-segmented positive-sense RNA virus.² It is spread by human-to-human transmission via droplets or direct contact. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. Among patients with pneumonia caused by SARS-CoV-2, fever was the most common symptom, followed by cough.³ Influenza, commonly called “the flu”, is an infectious disease caused by influenza viruses.⁴ The symptoms range from mild to severe and often include fever, runny nose, sore throat, muscle pain, headache, coughing, and fatigue. These symptoms typically begin 1-4 days after exposure to the virus and last for about 2-8 days.⁵

Test Principle

Panbio™ COVID-19/Flu A&B Rapid Panel (Nasopharyngeal) contains a membrane strip, which is pre-coated with immobilized anti-SARS-CoV-2, anti-Influenza A, and anti-Influenza B antibodies on the test line and mouse monoclonal anti-chicken IgY on the control line. Four types of gold conjugates (monoclonal antibodies specific to SARS-CoV-2, Influenza A, and Influenza B gold conjugates and chicken IgY gold conjugate) move upward on the membrane chromatographically and react with anti-SARS-CoV-2, anti-Influenza A, anti-Influenza B antibodies and pre-coated mouse monoclonal anti-chicken IgY antibody respectively. For a positive result, each pathogen of SARS-CoV-2, Influenza A and Influenza B will form a test line in the result window. Neither the test line nor the control line are visible in the result window prior to applying the patient specimen. A visible control line is required to indicate a test result is valid.

Intended Use

The Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is an *in vitro* diagnostic rapid test for the qualitative detection of SARS-CoV-2 and Influenza A&B antigen (Ag) in human nasopharyngeal swab specimens from individuals, who meet COVID-19 and/or Influenza clinical and/or epidemiological criteria.

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is for professional use and is intended to be used as an aid in the diagnosis of SARS-CoV-2 and Influenza infection. The product may be used in laboratory and non-laboratory environment.

The test provides preliminary test results. Negative results don't preclude SARS-CoV-2 and Influenza A&B infection and they cannot be used as the sole basis for treatment or other management decisions. Negative results should be combined with clinical observations, patient history, and epidemiological information.

Materials Provided

- 10 Test devices with desiccant in individual foil pouches
- 10 Extraction buffer tubes in a foil pouch
- 10 Nozzle caps
- 10 Sterilized nasopharyngeal swabs for specimen collection
- 1 Tube rack
- 1 Instructions for Use with Quick Reference Guide

Materials Required but not Provided

- Personal Protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves)
- Timer
- Biohazard container

Active Ingredients of Main Components

- Test device
Gold conjugate: monoclonal antibodies to SARS-CoV-2 Ag, anti-Influenza A Ag, and B Ag - gold colloid and Chicken IgY - gold colloid, Test line: Mouse monoclonal anti-SARS-CoV-2, Mouse monoclonal anti-Influenza A and Fragment antigen-binding anti-Influenza B, Control line: Mouse monoclonal anti-Chicken IgY
- Buffer: Tricine, Sodium Chloride, Tween 20, Tergitol 15-S-9, Sodium Azide (<0.1%)

Storage and Stability of Test Kit

1. The test kit should be stored at a temperature between 2 - 30 °C. Do not freeze the kit or its components.
2. When stored in a refrigerator, all kit components must be brought to room temperature (15 - 30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.
3. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use the test kit beyond its expiration date.
5. The shelf life of the kit is indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. Direct swab specimens should be tested immediately after collection. If immediate testing is not possible, the swab specimen can be kept in an extraction buffer tube filled with extraction buffer at room temperature (15 - 30°C) for up to 2 hours prior to testing.

Warnings

1. For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
2. These instructions must be strictly followed by a trained healthcare professional to achieve accurate results. All users have to read the instructions prior to performing a test.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen and buffer.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. swab, extraction buffer tube, test device) in a biohazard container as if they were infectious waste and dispose according to applicable local regulations.
8. Do not mix or interchange different specimens.
9. Do not mix reagent of those for other product.
10. Do not store the test kit in direct sunlight.
11. To avoid contamination, do not touch the head of provided swab when

- opening the swab pouch.
- The sterilized swabs should be used only for nasopharyngeal specimen collection.
 - To avoid cross-contamination, do not reuse the swabs for specimen collection.
 - Do not dilute the collected swab with any solution except for the provided buffer.
 - The buffer contains <0.1% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with a large volume of water.⁶

Test Procedure (Refer to Figures)

Nasopharyngeal Swab Specimens

Note: Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.

Test Preparation

- When stored in a refrigerator, all kit components must be brought to room temperature (15 - 30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.
- Remove the test device from the foil pouch. Place on a flat, horizontal and clean surface.
- Keep extraction buffer tube upright and remove foil. Place the extraction buffer tube in the tube rack.

Nasopharyngeal Specimen Collection & Extraction

- Tilt the patient's head back slightly about 45°-70° to straighten the passage from the front of the nose.
- Insert the swab through the nostril parallel to the palate.
 **Caution:** Use dedicated nasopharyngeal swab for specimen collection, provided in the kit.
- Swab should reach depth equal to distance from nostrils to outer opening of the ear.
 **Caution:** If resistance is encountered during insertion of the swab, remove it and attempt insertion in the opposite nostril.
- Gently rub and roll the swab, 5 times. Leave the swab in place for several

seconds to absorb secretions.

5. Slowly remove swab while rotating it and insert into the extraction buffer tube.
6. Swirl the swab tip in the buffer fluid inside the extraction buffer tube, pushing into the wall of the extraction buffer tube at least five times and then squeeze out the swab by squeezing the extraction buffer tube with your fingers.
7. Dispose of the used swab in accordance with your local regulations and biohazard waste disposal protocol.
8. Press the nozzle cap tightly onto the extraction buffer tube.

Reaction with Test Device

1. Dispense 4 drops of extracted specimen vertically into the specimen well on the device.
 **Caution:** Do not handle or move the test device until the test is complete and ready for reading.
2. Start timer. Read result at 15 minutes. Do not read result before 15 minutes or after 20 minutes.
3. Close the nozzle using cap and dispose of the extraction buffer tube according to your local regulations and biohazard waste disposal protocol.
4. Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

Test Limitations

1. The contents of this kit are to be used for the professional and qualitative detection of SARS-CoV-2 and Influenza A&B antigens from nasopharyngeal swab. Other specimen types may lead to incorrect results and must not be used.
2. Failure to follow the instructions for test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
3. A negative test result may occur if the specimen was collected, extracted or transported improperly. A negative test result does not eliminate the possibility of SARS-CoV-2 and Influenza A&B infection and should be confirmed by viral culture or a molecular assay.

4. Positive test results do not rule out co-infections with other pathogens.
5. Test results must be evaluated in conjunction with other clinical data available to the physician.
6. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.
7. Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is not intended to detect from defective (non-infectious) virus during the later stages of viral shedding that might be detected by PCR molecular tests.⁷
8. Positive results may occur in cases of infection with SARS-CoV or Albumin.
9. Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was not evaluated for the following influenza strains as referenced by WHO guideline: A/Darwin/9/2021 (H3N2), B/Austria/1359417/2021.
10. Wait 4 hours before repeating the test following an invalid result.
11. Clinical performance was evaluated with frozen samples, and test performance may be different with fresh samples.

Quality Control

Procedural Controls

The test device has a control line (“C”) and test lines (“FLU B”, “FLU A”, “COV 19”) on the surface of the test device. Neither the test line nor the control line are visible in the result window before applying a specimen. The control line is used for procedural control and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

External Controls

Good laboratory practice recommends the use of external control materials. Users should follow federal, state and local guidelines concerning the frequency of external control testing. Specific positive controls can be ordered by the following catalogue number: 62FK11C.

Performance Characteristics

1. Clinical performance with NP Specimen of Panbio™ COVID-19/Flu A&B Rapid Panel (NP)

The Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is designed for nasopharyngeal specimen collection from one nostril. The performance of Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was established with Nasopharyngeal swabs tested in Korea in September 2021 with 100 COVID-19 positives, 50 Influenza A positives, 50 Influenza B positives, and 200 negative specimens for COVID-19 and Influenza A&B antigen.

Samples were collected from subjects who presented with symptoms of Influenza A&B and/or COVID-19. One specimen was collected for testing on Panbio™ COVID-19/Flu A&B Rapid Panel (NP) and another for PCR. Swabs were collected and extracted into viral transport media. Samples were frozen after collection and thawed before testing. The performance of the Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was evaluated by comparison to a CE-marked RT-PCR method.

The study showed a sensitivity of 88.0% (95% CI: 80.0-93.6%) and specificity of 100.0% (95% CI: 98.2-100.0%) for SARS-CoV-2, a sensitivity of 92.0% (95% CI: 80.8-97.8%) and specificity of 100.0% (95% CI: 98.2-100.0%) for Influenza A, and a sensitivity of 100.0% (95% CI: 92.9-100.0%) and specificity of 100.0% (95% CI: 98.2-100.0%) for Influenza B.

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) Results for NP COVID-19

SARS-CoV-2		RT-PCR			Overall Sensitivity (95% CI)	Overall Specificity (95% CI)
		Positive	Negative	Total		
Panbio™ COVID-19/Flu A&B Rapid Panel (NP)	Positive	88	0	88	88.0% (80.0-93.6%)	100.0% (98.2-100.0%)
	Negative	12	200	212		
	Total	100	200	300		

- Positive agreement of the Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is higher with specimens of Ct values ≤ 30 with a sensitivity of 94.4% (95%CI:87.5-98.3%) and lower with specimens of Ct values >30 with a sensitivity of 30.0% (95%CI:6.7-65.3%) for SARS-CoV-2. As suggested in references 7-9, patients with Ct values >30 are no longer contagious.^{7,8,9}

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) Results for NP Influenza A

Influenza A		RT-PCR			Overall Sensitivity (95% CI)	Overall Specificity (95% CI)
		Positive	Negative	Total		
Panbio™ COVID-19/Flu A&B Rapid Panel (NP)	Positive	46	0	46	92.0% (80.8-97.8%)	100.0% (98.2-100.0%)
	Negative	4	200	204		
	Total	50	200	250		

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) Results for NP Influenza B

Influenza B		RT-PCR			Overall Sensitivity (95% CI)	Overall Specificity (95% CI)
		Positive	Negative	Total		
Panbio™ COVID-19/Flu A&B Rapid Panel (NP)	Positive	50	0	50	100.0% (92.9-100.0%)	100.0% (98.2-100.0%)
	Negative	0	200	200		
	Total	50	200	250		

2. Detection Limit

The detection limit of the Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was determined on the following virus strains:

- 1.54×10^2 TCID₅₀/ml of SARS-CoV-2 which was isolated from a COVID-19 confirmed patient in Korea.
- 4.07×10^2 TCID₅₀/ml of Influenza A H1N1/California/07/09 strain.
- 1.14×10^2 TCID₅₀/ml of Influenza A H3N2/Perth/16/09.
- 2.46×10^4 IU/ml of Influenza B/Wisconsin/1/2010.

This influenza A viral strains were obtained from ZeptoMetrix with titer information. The performance characteristics for influenza A virus can vary depending on the subtypes, region, and year of occurrence since influenza A viruses have undergone minor amino acid changes in the target epitope region.

3. Hook Effect

There is no hook effect at:

- $1.0 \times 10^{5.8}$ TCID₅₀/ml with SARS-CoV-2 culture fluid which was isolated from a COVID-19 confirmed patient in Korea.
- 1.17×10^5 TCID₅₀/ml with Influenza A H1N1/California/07/09.
- 1.17×10^5 TCID₅₀/ml with Influenza A H3N2/Perth/16/09.
- 1.40×10^7 IU/ml with Influenza B/Wisconsin/1/2010.
- $5.0 \times 10^{7.25}$ TCID₅₀/ml with Influenza A H1N1/New Jersey/8/76.
- 3.40×10^9 TCID₅₀/ml with Influenza A/Denver/1/57.
- $5.0 \times 10^{6.5}$ TCID₅₀/ml with Influenza A/Hong Kong/8/68.

4. Cross Reactivity

Cross-reactivity of Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was evaluated by testing 76 viruses and 37 other microorganisms (5 replicates of each). The human Influenza A strains, as expected, reacted with the Influenza A test line and did not cross react with the influenza B or COVID-19 test lines. The human Influenza B strains, as expected, reacted with the Influenza B test line and did not cross react with the influenza A or COVID-19 test lines. The non-human influenza strains did not react with the Influenza A, B or COVID-19 test lines. Panbio™ COVID-19/Flu A&B Rapid Panel (NP) has cross-reactivity with Human-SARS-coronavirus Nucleoprotein at a concentration of ≥ 200 pg/ml due to SARS-CoV having a high homology to the SARS-CoV-2 (79.6%).

Tested Viruses		
SARS-CoV Nucleoprotein	Influenza A H1N1 (California/07/09)	Influenza A H1N1 (A/Beijing/262/1995)
Influenza A H1N1 (A/New Jersey/8/76)	Influenza A H1N1 (A/Denver/1/57)	Influenza H1N1 (A/Puerto Rico/8/34)
Influenza A H1N1 (A/ FM/1/47)	Influenza H3N2 (A/Hong Kong/8/68)	Influenza A H3N2 (A/Victoria/3/75)
Influenza A H3N2 (A/Port Chalmers/1/73)	Influenza A H3N2 (Perth/16/09)	Influenza A H3N2 (Wisconsin/67/05)
Influenza A H5N1 (A/ China/2006)	Influenza A H5N2 (A/chicken/ Iowa/ 04-20/2015)*	Influenza A H7N1 (A/turkey/ Italy/4602/99)*
Influenza A H7N2 (A/ ruddy turnstone/ New Jersey/563/2006)*	Influenza A H7N7 (A/chicken/ Netherlands/1/03)*	Influenza B (B/Flori- da/78/2015, Victoria Lineage)
Influenza B (B/ Florida/4/2006, Yamagata Lineage)	Influenza B (B/Maryland/1/59, Victoria Lineage)	Influenza B (B/GL/1739/54)
Influenza B (B /Taiwan/2/62)	Influenza B (B/Malaysia/2506/04, Victoria Lineage)	Influenza B (B/Wisconsin/1/2010, Yamagata Lineage)
Influenza A H1N1 (A/Wisconsin/588/2019)/ (A/Victoria/2570/2019)	Influenza A H3N2 (A/Cambodia/ e0826360/2020)	Influenza B (B/ Phuket/3073/2013, Yamagata Lineage)
Influenza B (B/ Washington/02/2019, Victoria Lineage)		

* non-human influenza strain

There was no cross-reaction with the viruses and microorganisms listed below.

Tested Viruses		
HCoV-HKU1	HCoV-NL63	HCoV-OC43
MERS-CoV	HCoV-229E	Human adenovirus 1
Human adenovirus 2	Human adenovirus 3	Human adenovirus 4
Human adenovirus 5	Adenovirus 6	Adenovirus 7
Adenovirus 10	Human adenovirus 18	Adenovirus 21
Norovirus	Epstein Barr Virus	Human herpesvirus (HSV) 1
Human herpesvirus (HSV) 2	Human cytomegalovirus	Mumps virus
Parainfluenza virus type 4A	Parainfluenza virus type 4B	Parainfluenza type 1
Parainfluenza type 2	Parainfluenza type 3	Human respiratory syncytial virus A
Human respiratory syncytial virus B	Measles virus	hMPV 16 type A1
Human rhinovirus type 1B	Rhinovirus type A2	Human rhinovirus 14
Rhinovirus A16	Rotavirus A	Sendai virus
Rubella virus	Coxsackievirus A9	Coxsackievirus B4
Coxsackievirus B5	Human echovirus 2	Human echovirus 6
Echovirus 7	Echovirus 11	Human enterovirus 71
Poliovirus 1	Enterovirus D68	

Tested Microorganisms		
<i>Bordetella pertussis</i>	<i>Mycobacterium avium</i>	<i>Staphylococcus epidermidis</i>
<i>Candida albicans</i>	<i>Mycobacterium tuberculosis</i>	<i>Staphylococcus haemolyticus</i>
<i>Chlamydomphila pneumoniae</i>	<i>Mycoplasma pneumoniae</i>	<i>Staphylococcus saprophyticus</i>
<i>Enterococcus faecalis</i>	<i>Neisseria lactamica</i>	<i>Streptococcus pneumoniae</i>
<i>Escherichia coli</i>	<i>Neisseria gonorrhoeae</i>	<i>Streptococcus pyogenes</i>
<i>Fusobacterium necrophorum</i>	<i>Neisseria meningitidis</i>	<i>Streptococcus salivarius</i>
<i>Haemophilus influenzae</i>	<i>Neisseria sp.</i> (<i>Neisseria lactamica</i>)	<i>Streptococcus anginosus</i>
<i>Klebsiella pneumoniae</i>	<i>Proteus vulgaris</i>	<i>Pneumocystis jirovecii</i> (PJP)
<i>Legionella pneumophila</i>	<i>Pseudomonas aeruginosa</i>	Pooled human nasal wash
<i>Moraxella catarrhalis</i>	<i>Corynebacterium diphtheriae</i>	<i>Enterococcus faecium</i>
<i>Lactobacillus plantarum</i>	<i>Neisseria sicca</i>	<i>Streptococcus agalactiae</i>
<i>Streptococcus bovis</i>	<i>Streptococcus oralis</i>	<i>Arcanobacterium haemolyticum</i>
<i>Staphylococcus aureus</i>		

5. Interfering Substances

The following 53 potentially interfering substances (5 replicates of each) have no impact on Panbio™ COVID-19/Flu A&B Rapid Panel (NP) with the exception of Albumin at concentrations >55 g/L.

Type	Tested Substances			
Endogenous	Mucin	Human blood	Bilirubin	Antinuclear antibody (ANA)
	Blood lipid	Hemoglobin	Rheumatoid factor	Human Anti-Mouse Antibody
	Albumin	Mometasone	Lopinavir	Acetylsalicylic acid
Exogenous	Phenylephrine	Fluticasone	Ritonavir	Ibuprofen
	Oxymetazoline	Chlorpheniramine	Arbidol	Sulphur (Sulfur)
	Sodium chloride	Diphenhydramine	Levofloxacin	Luffa operculata
	Nasal spray (NasalCrom)	Histamine dihydrochloride	Azithromycin	Benzocaine
	Nasal spray (Equate)	Histamine	Meropenem	Menthol
	Beclomethasone	Galphimia glauca	Tobramycin	Acetaminophen
	Dexamethasone	Histaminum hydrochloricum	Mupirocin	Mouthwash (Crest)
	Flunisolide	Zanamivir	Albuterol	Mouthwash (Listerine)
	Triamcinolone acetoneide	Ribavirin (Rebetol)	Ephedrine	Mouthwash (Therabreath)
	Budesonide	Tamiflu (Oseltamivir)	Guaicol glyceryl ether	Heparin
	Phenylpropanolamine			

6. Repeatability & Reproducibility

Repeatability & Reproducibility of Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was established using 3 positive specimens (moderate, weak, LoD proximity concentration) per each strain (SARS-CoV-2, Influenza A H1N1, Influenza A H3N2, and Influenza B) and negative specimen. There were no differences observed within-run, between-run, between-lots, between-sites, and between-days.

PREPARATION

- 1 When stored in a refrigerator, all kit components must be brought to room temperature (15 - 30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.

Note: Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.

2 Open the package and look for the following:

1. Test device with desiccant in individual foil pouch
2. Extraction buffer tube in the foil pouch
3. Nozzle cap
4. Sterilized nasopharyngeal swab for specimen collection
5. Tube rack
6. Instructions for Use with Quick Reference Guide


3 Carefully read these instructions prior to using Panbio™ COVID-19/Flu A&B Rapid Panel (NP).

4 Look at the expiration date of the kit box. If the expiration date has passed, select another kit for use.

5 Open the foil pouch of test device and look for the following:


1. Result window
2. Specimen well

Then, label the device with the patient identifier.


TEST PROCEDURE


- 1 Keep extraction buffer tube upright and remove foil. Place the extraction buffer tube in the tube rack.

- 2 Tilt the patient's head back (about 45°-70°). Insert the swab through the nostril. Gently rub and roll the swab, 5 times. Leave the swab in place for several seconds. Slowly remove swab.

- 3 Insert the swab specimen into the extraction buffer tube. Swirl the swab tip in the buffer fluid inside the extraction buffer tube, pushing into the wall of the extraction buffer tube at least five times and then squeeze out the swab by squeezing the extraction buffer tube with your fingers.

TEST PROCEDURE

- 4** Dispose of the used swab in accordance with your local regulations and biohazard waste disposal protocol.


- 5** Press the nozzle cap tightly onto the extraction buffer tube.

- 6** Dispense 4 drops of extracted specimen vertically into the specimen well on the device.

 **Caution:** Do not handle or move the test device until the test is complete and ready for reading.

TEST PROCEDURE

7 Start timer. Read result at 15 minutes.

Do not read result after 20 minutes.


Note: A Control (“C”) line may appear in the result window in a few minutes but a test line (“FLU B”, “FLU A”, “COV 19”) may take as long as 15 minutes to appear.

WAIT 15 MINUTES

READ AT 15-20 MINUTES





8 Close the nozzle using cap and dispose of the extraction buffer tube according to your local regulations and biohazard waste disposal protocol.

TEST INTERPRETATION

POSITIVE

- **COVID-19 positive:** The presence of two colored lines (test line “COV 19” and control line “C”) within the result window.
- **Influenza A positive:** The presence of two colored lines (test line “FLU A” and control line “C”) within the result window.
- **Influenza B positive:** The presence of two colored lines (test line “FLU B” and control line “C”) within the result window.
- **Multi positives:** The presence of three or four colored lines (test line “FLU B”, “FLU A”, “COV 19” and control line “C”) within the result window.

 **Caution:** The presence of any test line (“FLU B”, “FLU A”, “COV 19”), no matter how faint, indicates a positive result.


Note: In the case of multi-positives, User should be aware of this unique situation as both can present with similar symptoms but vary in treatment.

Note that the test provides preliminary test results.

TEST INTERPRETATION

NEGATIVE


The presence of only the control line (“C”) and no test line (“FLU B”, “FLU A”, “COV 19”) within the result window indicates a negative result.

INVALID

If the control line (“C”) is not visible within the result window after performing the test, the result is considered invalid (even though the test line may be visible).

Instructions may not have been followed correctly. It is recommended to wait 4 hours, read the IFU again and collect a new specimen and re-test using a new test device.

- 9 Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

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1. Rothan HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. *J Autoimmun.* 2020 May;109:102433. doi: 10.1016/j.jaut.2020.102433.
2. Guo YR, Cao QD, Hong ZS, et al. The origin, transmission and clinical therapies on coronavirus disease 2019 (COVID-19) outbreak - an update on the status. *Mil Med Res.* 2020 Mar 13;7(1):11. doi: 10.1186/s40779-020-00240-0.
3. Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. *J Antimicrob Agents.* 2020 Mar;55(3):105924. doi: 10.1016/j.jantimicag.2020.105924. Epub 2020 Feb 17.
4. "Flu Symptoms & Diagnosis". Centers for Disease Control and Prevention (CDC). 10 July 2019. Retrieved 24 January 2020.
5. "Flu Symptoms & Complications". Centers for Disease Control and Prevention (CDC). 26 February 2019. Retrieved 6 July 2019.
6. Current Intelligence Bulletin 13: Explosive Azide Hazard DHHS (NIOSH) Publication Number 78-127 August 16, 1976.
7. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
8. CDC. Duration of Isolation of Precautions for Adults with COVID-19. (2020).
9. Bullard, et al. Predicting Infectious Severe Acute Respiratory Syndrome Coronavirus 2 From Diagnostic Samples. *Clin Infect Dis.* 2020 Dec 17;71(10):2663-2666. doi: 10.1093/cid/cia638.

GLOSSARY OF SYMBOLS / GLOSARIO DE SÍMBOLOS / GLOSSÁRIO DE SÍMBOLOS / GLOSAR DE SIMBOLURI / REČNIK SIMBOLA / УСЛОВНЫЕ ОБОЗНАЧЕНИЯ / OBJAŚNIENIE SYMBOLI / SLOVNÍK SYMBOLŮ

	Temperature limitation Limitación de temperatura Limitação de temperatura Limite de temperatură Ograničenje temperature Ограничение температуры Ograniczenie temperatury Omezení týkající se teploty
	For <i>in vitro</i> diagnostic use only Para uso diagnóstico <i>in vitro</i> únicamente Somente para uso em diagnóstico <i>in vitro</i> Destinat utilizării doar pentru diagnostic <i>in vitro</i> Samo za <i>in vitro</i> dijagnostičku upotrebu Только для диагностики <i>in vitro</i> Wyłącznie do użytku w diagnostyce <i>in vitro</i> Pouze k použití pro diagnostiku <i>in vitro</i>
	Do not reuse No reutilizar Não reutilizar A nu se reutiliza Ne koristiti ponovo Повторное использование запрещено Nie używać ponownie Nepoužívejte opakovaně

	Do not use if package is damaged No lo use si el envase está dañado Não usar se a embalagem estiver danificada A nu se utiliza dacă ambalajul este deteriorat Ne koristiti ako je pakovanje oštećeno Не использовать при повреждении упаковки Nie używać, jeśli opakowanie jest uszkodzone Nepoužívejte, je-li je obal poškozen
	Lot Number Número de lote Número do lote Număr de lot Broj serije Номер партии Numer partii Číslo šarže
	Catalog Number Número de catálogo Número de catálogo Număr de catalog Kataloški broj Номер по каталогу Numer katalogowy Katalogové číslo
	Consult instructions for use Consulte las instrucciones de uso Consulte as instruções de uso A se consulta instrucțiunile de utilizare Pogledajte uputstvo za upotrebu Обратитесь к инструкциям по применению Sprawdzić w instrukcji użycia Prostudujte si návod k použití

	Keep dry Manténgase seco Mantenha seco A se păstra uscat Održavajte suvim Беречь от влаги Chronić przed wilgocią Uchovávejte v suchu
	Biological Risks Riesgos biológicos Riscos biológicos Riscuri biologice Biolóški rizici Биологическая опасность Zagrożenie biologiczne Biologická rizika
	Use By Fecha de caducidad Usar até A se utiliza până la Upotrebiti do Срок годности Użyć przed Spotřebujte do
	Manufacturer Fabricante Fabricante Producător Proizvodač Производитель Producent Výrobce

	Date of manufacture Fecha de fabricación Data de fabrica��o Data fabrica��iei Datum proizvodnje ��ata ��izgotovljeniya Data produkcji Datum v��roby
	Keep away from sunlight Mant��ngase alejado de la luz solar Mantenha longe da luz solar A se feri de lumina solar�� Dr��zati dalje od sun��eve svetlosti ��ere��y ot soln��chnykh luchey Chroni�� przed ��wiatfem s��lonecznym Udr��ujte mimo dosah slune��n��ho z��ření
	CE mark Marca CE Marca CE Marcaj CE CE oznaka ��ar��irovka CE Znak CE Ozna��en�� CE
	Contains sufficient for X tests Contiene suficiente producto para X pruebas Cont��m o suficiente para X testes Con��ine o cantitate suficient�� pentru X teste Sadr��aj dovoljno za X testova ��oder��imogo ��ostato��no ��lya provedeniya X testov Zawiera ilo��c wystarczaj��c�� do wykonania X test��w Obsahuje dostate��n�� mno��stv�� pro X test��

	Caution Precaución Cuidado Atenție Oprez Внимание Ostrzeżenie Upozornění
	Sterilized using ethylene oxide Esterilizado con óxido de etileno Esterilizado com óxido de etileno Sterilizat cu oxid de etilenă Sterilisano etilen-oksidad Стерилизовано с помощью этиленоксида Wyrób sterylizowany tlenkiem etylenu Sterilizováno etylénoxidem
	Sterilized using irradiation Esterilizado mediante irradiación Esterilizado por irradiação Sterilizat prin iradiere Sterilisano zračenjem Стерилизовано путем облучения Wyrób sterylizowany za pomocą promieniowania Sterilizováno ozařováním
	Do not re-sterilize No lo vuelva a esterilizar Não reesterilizar A nu se resteriliza Ne sterilisati ponovo Не стерилизовать повторно Nie sterylizować ponownie Nesterilizujte opakovaně

TECHNICAL SUPPORT / SOPORTE TÉCNICO / SUPORTE TÉCNICO / ASISTENȚĂ TEHNICĂ / TEHNIČKA PODRŠKA / ТЕХНИЧЕСКАЯ ПОДДЕРЖКА / WSPARCIE TECHNICZNE / TECHNICKÁ PODPORA	
Europe and Middle East / Europa y Medio Oriente / Europa e Oriente Médio / Europa și Orientul Mijlociu / Evropa i Bliski istok / Европа и Ближний Восток / Europa i Bliski Wschód / Evropa a Střední východ	+44 161 483 9032 EME.TechSupport@abbott.com
Africa / África / África / Africa / Afrika / Африка / Afryka / Afrika	+27 10 500 9700 arcis.techsupport@abbott.com
Russia and CIS / Rusia y CEI / Rússia e CEI / Rusia și CSI / Rusija i ZND / Россия и СНГ / Rosja i WNP / Rusko a SNS	+7 499 403 9512 arcis.techsupport@abbott.com
Asia Pacific / Asia Pacífico / Ásia-Pacífico / Asia Pacific / Azija-Pacifik / Азиатско-Тихоокеанский регион / Azja i rejon Pacyfiku / Asie a Tichomoří	+61 7 3363 7711 AP.TechSupport@abbott.com
Latin America / América Latina / América Latina / America Latină / Latinska Amerika / Латинская Америка / Ameryka Łacińska / Latinská Amerika	+57 60 1794 5968 LA.TechSupport@abbott.com


Canada / Canadá / Canadá / Canada / Kanada / Канада / Kanada / Kanada	+1 800 818 8335 Canproductsupport@abbott.com
---	---

- MEMO -

- MEMO -

- MEMO -

- MEMO -

Product - COVID-19/Flu (Cluster A) / IFU Size - 110 x 140mm	 Black 100	Rev. No. 62FK10A-08-A0 Rev. Date - 2022. 04. 26
--	--	---