3. Study Site 3
259 patients serum specimens collected from Sri Lanka were tested on the SD BIOLINE Dengue Duo at the Mahidol-Oxford Tropical Medicine Unit (MORU), Bangkok, Thailand. The specimens were comprised of 99 acute dengue specimens and 160 non-dengue specimens with other confirmed febrile illness. 99 acute dengue specimens were confirmed by RT-PCR and MAC/GAC ELISA. The results were used to determine the sensitivity and specificity.
Result (% CI): Sensitivity: 92.9% (83.9 - 97.1%), Specificity: 88.8% (82.8 - 93.2%).

6. Cross-reactivity test with other flavivirus mediated and mosquito-borne disease
The following 4 potential cross-reacting pathogens had no effect on test results of SD BIOLINE Dengue Duo. Even though there are no cross reactivity in the results of internal study, the possibility can't be excluded completely.
Disease **Dengue/IgG/Positive/Total** **Dengue/IgG/Positive/Total** **Dengue/IgA/Positive/Total**

Japanese Encephalitis	0/25	0/25	0/25
Yellow Fever	0/25	0/25	0/25
Malaria P. falciparum	0/25	0/25	0/25
Malaria P. vivax	0/25	0/25	0/25
Total	0/100	0/100	0/100

7. Reproducibility of SD BIOLINE Dengue Duo test results has been demonstrated by within-run, between-run and lot-to-lot studies using in-house reference panels. All values were identical to reference panel across all criteria.

8. The performance of SD BIOLINE Dengue Duo kit was evaluated with potential relevant interfering substances. Hemolytic, lipemic, icteric specimens and those containing rheumatoid factors were investigated. None of the specimens were found to interfere with this test kit.

Product Disclaimer :
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The safety of the diagnosis should consult a doctor for further confirmation of the test result.

Warning :
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect positive or negative diagnosis using this product.

4. Study Site 4
Sensitivity of SD BIOLINE Dengue Duo according to the type and phase of dengue infection
A total of 397 serum specimens collected from Mexico were tested on the SD BIOLINE Dengue Duo in order to evaluate the performance. The specimens were confirmed as having dengue infection by ELISA test. On the basis of NS IgG and IgA ELISA results and clinical date of patient, 310 specimens were diagnosed with dengue virus infection and 87 specimens were negative for all 3 ELISA. The performance of SD BIOLINE Dengue Duo was evaluated against the overall results of ELISA test.

Specimen type	Specimen	SD BIOLINE Dengue Duo Result (%) CI
Primary acute dengue	122	110 (86.8 - 95.8%)
Primary convalescent dengue	17	11 (64.71% (39.05 - 59.36%))
Secondary acute dengue	122	117 (59.50% (91.97 - 99.83%))
Secondary convalescent dengue	43	43 (87.75% (77.55 - 97.95%))
Non dengue	87	9 (78) Specificity: 89.6% (82.68 - 96.63%)

* Dengue positive if at least 1 of NS1 IgG or IgM marker was positive
* Non dengue if all 3 markers were negative

5. Study Site 5
A total of 199 serum specimens collected and tested on the SD BIOLINE Dengue Duo at Ministerio de Salud (MINSa), Argentina. The serum specimens were comprised of 149 positive and 50 negative specimens confirmed by ELISA and RT-PCR.

Specimen type	Days after symptom onset	No. of specimens	SD BIOLINE Dengue Duo
Positive	<3	49	24 Positive 25 Negative
	4-6	50	35 Positive 15 Negative
	>6	50	50 Positive 0 Negative
Negative	50	0	50
Sensitivity (95% CI)			73.2% (65.5 - 79.6%)
Specificity (95% CI)			100% (92.9 - 100%)

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