# DETERMINE™ HIV EARLY DETECT

Abbott

**REF** 7D2842, 7D2843

January 2020  
241851/R9

1

**Key to symbols used/**  
**Clave de los símbolos utilizados/**  
**Légende des symboles utilisés/**  
**Designação dos símbolos utilizados**

**30°C** Store at 2-30°C/  
Guardar a temperaturas entre 2 y 30°C/  
Conserver entre 2-30°C/  
Armazenar a 2-30°C

**REF**

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## Advice Line

For further information, please contact your distributor, or call to one of the following Abbott Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+ (44) 161 483 9032	EMEproductsupport@abbott.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@abbott.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCIProductsupport@abbott.com
Latin America	+ (57) 2 661 8797	LAproductsupport@abbott.com

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This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

## NAME AND INTENDED USE

**Determine™ HIV Early Detect** is an *in vitro*, visually read, qualitative immunoassay for the detection of antibodies (Ab) to HIV-1 and HIV-2 and the detection of non-immunocomplexed (free) HIV-1 p24 antigen (Ag) in human capillary and venous whole blood, plasma or serum. Therefore, any reactivity on 1) the Ab bar alone 2) the Ag bar alone or 3) both the Ab and Ag bars simultaneously, is considered a reactive result suggestive of infection with HIV. The test is for professional use only and intended to be used as an aid in the diagnosis of HIV-1/HIV-2 infection in adults and children over the age of 18 months.

## SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of CD4 helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The multiplication of the HIV in the infected cells leads to cell rupture and thus the release of HIV virus particles, which are first detected in the form of HIV RNA and next in the form of HIV antigen.<sup>1,2</sup> This is followed by production of specific antibodies to either HIV-1 or HIV-2.<sup>3,4,5</sup> HIV antigen may become undetectable at this time because of the formation of antibody-antigen complexes.

## BIOLOGICAL PRINCIPLES OF THE PROCEDURE

**Determine HIV Early Detect** is an immunochromatographic test for the qualitative detection of free HIV-1 p24 antigen and antibodies to HIV-1 and HIV-2.

Specimen is added to the sample pad. The specimen mixes with biotinylated anti-p24 antibodies and selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens, synthetic HIV-2 peptide and anti p24 mouse monoclonal antibody. This mixture continues to migrate through the solid phase to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides at the Antibody (Ab) window, immobilized avidin at the Antigen (Ag) window. If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens and synthetic HIV-2 peptide and to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides, forming one red bar at the Ab window site. If antibodies to HIV-1 and HIV-2 are absent, the selenium colloid-conjugates flow past the Ab window and no red bar is formed at the Ab window site.

If free HIV-1 p24 antigen is present in the specimen, the antigen binds to the biotinylated anti-p24 antibodies and the selenium colloid-conjugate coated with anti p24 mouse monoclonal antibody. This complex binds to the immobilized avidin forming a red bar at the Ag window site. If HIV-1 p24 antigen is not present, both the biotinylated anti-p24 antibodies and selenium colloid-conjugate flow past the Ag window and no red bar is formed at the Ag window site.

To ensure assay validity, a procedural control bar is incorporated in the assay device at the Control window.

## CONTENTS

**Determine HIV Early Detect 20 Test (7D2842) or 100 Test (7D2843)**

- **Determine HIV Early Detect Test Card, 2 or 10 cards (containing 10 tests/card) coated with HIV-1/2 recombinant antigen, synthetic peptides, anti-p24 antibodies and avidin.**

## Accessories (required but not provided)

### For testing Whole Blood samples

- Chase Buffer (7D2243) 1 Bottle (2.5 mL) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

### For testing Whole Blood samples (fingerstick assay)

- EDTA Capillary Tubes (7D2222)

### **Determine HIV Early Detect SET (7D2843SET) (100 Test for testing whole blood samples)**

- **Determine HIV Early Detect (7D2843)**
- Chase Buffer (7D2243)
- EDTA Capillary Tubes (7D2222)
- Blood Lancet (7D2232 or 7D2233)

## Materials Required But Not Provided

- Disposable gloves, timing device
- Micropipette capable of delivering 50 µL (other than fingerstick)
- Alcohol swab, gauze pad, lancet (for fingerstick)

## WARNINGS AND PRECAUTIONS

### For *In Vitro* Diagnostic Use.

For professional use only.

Safety data sheet available for professional user on request.

**CAUTION:** When handling specimens and reagents, use appropriate biosafety practices<sup>7,8</sup>. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using suitable disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.<sup>7,8</sup>
- Decontaminate and dispose of all specimens, used test strips, and other potentially contaminated materials in accordance with local regulations.<sup>7,8</sup>

## STORAGE

Store **Determine HIV Early Detect** Test Cards and Chase Buffer at 2-30°C until expiration date.

- When handled and stored as directed, kit components are stable until the expiration date. Do not use kit components beyond expiration date.
- Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.
- Do not use devices that have become wet or if the packaging has become damaged.

## SPECIMEN COLLECTION

### Serum, Plasma, and Whole Blood Collection by Venipuncture.

Use EDTA collection tubes for whole blood and plasma specimens.

- Collect human serum, plasma and whole blood by venipuncture aseptically.
- To obtain serum, separate from the clot. To obtain plasma, separate from the packed cells. Separate specimens as soon as possible to avoid any hemolysis.

## Whole Blood Collection by Fingerstick<sup>6</sup> (See Fig.1)

Use EDTA Capillary Tubes (7D2222).

**CAUTION:** Glass capillaries may be damaged during transportation or when in use. Handle with care in order to avoid injury when removing from the package as well as during use and during disposal. Before collecting a fingerstick specimen, place a capillary tube on a clean dry surface.

1. Choose the fingertip of the middle, ring, or index finger (whichever is the least calloused). Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
2. Clean fingertip with alcohol; allow to air dry.
3. Position the hand palm-side up. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
4. Wipe away the first drop of blood with a sterile gauze pad.
5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the capillary tube to the drop of blood. Avoid air bubbles.

Fill the tube with blood between the 2 lines marked on the capillary tube (50 µL).

Fig.1

## SPECIMEN STORAGE

• Store serum and plasma specimens at 2-8°C and run the test within 7 days of collection. If testing is delayed more than 7 days, freeze the specimen at -20°C or colder.

- Avoid repeated freeze/thaw cycles.
- If serum or plasma specimens show particulate matter or turbidity, centrifuge at 10,000g for 5 minutes at room temperature before sampling. Carefully remove the 50 µL test sample from the supernatant. If a lipid layer is formed on the surface of the liquid, ensure to take the sample from the clear liquid below that layer.
- For whole blood collected by venipuncture, store at 2-8°C. Do not freeze whole blood specimens. Run the test within 7 days of collection.
- For whole blood collected by fingerstick, test immediately.

## TEST PROCEDURE

Remove the desired number of test strips from the 10-test card by bending and tearing at the perforation. The test is to be performed at 15-30°C.

### NOTE:

- To preserve the lot number which appears on the left side of the test card, remove individual test strips starting from the right side of the test card. The lot number and expiry date are not printed on the individual test strips.
- After removing the protective foil cover from each test strip, start the assay immediately. At under 80% or less relative humidity, the test can be used within 30 minutes after removing the protective foil cover.
- For serum or plasma, ensure thorough mixing of sample prior to use. For whole blood sample, mix well by gentle inversion of the tube immediately before testing.
- Running the test in high temperature/low humidity may affect the appearance of the Ag or Ab bar. If the test strip is partially dried and difficult to read at 20 to 40 minutes, the test should be repeated using a new test strip and result read at 20 minutes. When the test strip is partially dried, it appears as mixed white spot and grayish area.
- If serum or plasma sample does not flow or shows abnormal flow, such as stopping in the middle of the window, centrifuge the specimen and repeat the test with a new test strip.

**CAUTION:** Abnormal flow may occur with a whole blood (fingerstick) sample, if the capillary tube is not placed in the middle of the sample pad at an upright (vertical) position. If this occurs, collect a new sample and repeat the test using a new strip.

1. Remove the protective foil cover from each test.

### For serum or plasma samples:

- a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
- b. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.

### For whole blood (venipuncture) samples:

- a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).

b. Wait one minute to allow the sample to be absorbed, then apply one drop of Chase Buffer to the sample pad, holding the bottle vertically.

c. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.

### For whole blood (fingerstick) samples using an EDTA Capillary Tube:

a. Place the capillary tube containing the blood sample onto the middle of the sample pad (marked by the arrow symbol) at an upright (vertical) position.

b. Wait until all the blood is transferred from the capillary tube to the sample pad. Then immediately apply one drop of Chase Buffer to the sample pad, holding the bottle vertically.

**Caution:** do not lift the capillary tube from the sample pad before all the blood has been transferred – a bubble may form which will prevent the complete transfer of sample and invalidate the test. It may take more than one minute for full transfer of the sample.

c. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.

## QUALITY CONTROL

To ensure assay validity, a procedural control is incorporated in the device and is labeled "C". If the control bar does not appear by assay completion, the test result is invalid. Repeat the test using a new test strip.

## INTERPRETATION OF RESULTS (See pictures)


### NOTES:

- Any reactivity on 1) the Ab bar alone 2) the Ag bar alone or 3) both, the Ab and Ag bars simultaneously, is considered a reactive result suggestive of infection with HIV.
- Interpret any visible red bar (even very faint) in the control window as a valid result.
- The test result is reactive even if the AG or AB bars appear lighter or darker than the control bar. Any visible red bar, no matter how faint, is interpreted as reactive.
- A test which gives very high background should be considered invalid.
- If an invalid result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

## ANTIBODY REACTIVE

### (Two Bars - Control and AB Bars)

Red bars appear in both the control window (labeled "C") and in the Ab window (labeled "AB") of the strip. Interpret any visible red bar, no matter how faint, in the Ab window as reactive.

### (Antigen (p24) REACTIVE (Two Bars - Control and AG Bars))

Red bars appear in both the control window (labeled "C") and in the Ag window (labeled "AG") of the strip. Interpret any visible red bar, no matter how faint, in the Ag window as reactive. The presence of only an antigen response suggests that the infection is at an early stage. Follow up testing may be suggested in order to track the expected future detection of antibodies.

### ANTIBODY REACTIVE AND ANTIGEN (p24) REACTIVE (Three Bars - Control, AB and AG Bars)

Red bars appear in the control window (labeled "C"), the Ab window (labeled "AB") and the Ag window (labeled "AG") of the strip. Interpret any visible red bar in the Ab and Ag windows as reactive.

### NON-REACTIVE (One Bar - Control Bar)

One red bar appears in the control window of the strip (labeled "C"), and no red bar appears in the Ab and Ag windows of the strip (labeled "AG" and "AB").

### INVALID (No Control Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the Ab or Ag window of the strip, the result is invalid. Repeat the test using a new test strip.

## LIMITATIONS OF THE PROCEDURE

- **Determine HIV Early Detect** is designed to detect antibodies (Ab) to HIV-1 and/or HIV-2 and non-immunocomplexed (free) HIV-1 p24 antigen (Ag), in human serum, plasma and capillary and venous whole blood specimens. Other body fluids or pooled specimens may not give accurate results and should not be used.
- The intensity of the AB and AG bars does not necessarily correlate to the titer of antibody and antigen in the specimen, respectively.
- A reactive result for antibodies to HIV-1/2 combined with a non-reactive result for HIV-1 p24 antigen does not preclude the possibility of acute HIV infection.
- Reactivity for both antibody and antigen lines or for antigen reactivity only may be confirmed by 4<sup>th</sup> generation EIA or NAT where the IVD is indicated as an aid for diagnosis. Antibody reactivity only can also be confirmed by using another HIV Ab RDT. If non-reactive in subsequent testing, retest at least 14 days after the original test is recommended. Follow the local guideline for confirmatory methods if applicable.
- No test provides absolute assurance that a specimen does not contain low levels of HIV-1 p24 antigen and/or antibodies to HIV-1 and HIV-2 such as those present at a very early stage or late stage of HIV infection.
- A non-reactive result for both antibodies to HIV-1/2 and HIV-1 p24 antigen does not preclude the possibility of exposure to or infection with HIV-1 or HIV-2 viruses.
- The absence of AG bar may occur when all p24 antigen is bound by antibodies. When high levels of antibodies against the p24 antigen are present in the blood after seroconversion, the antibodies tend to bind to the antigens, forming immunocomplexes. **Determine HIV Early Detect** detects only non-immunocomplexed (free) antigens; it does not detect immunocomplexed (bound) antigens.
- Some known HIV-infected persons taking antiretroviral medication have been shown to produce false negative results when tested by rapid diagnostic tests.<sup>10,11,12</sup>
- Where clinical presentation is inconsistent with the non-reactive **Determine HIV Early Detect** test result, the individual should be re-tested for antibodies to HIV with **Determine HIV Early Detect** at least 14 days after the original test, or with a NAT assay that is indicated for diagnosis.
- Whole blood or plasma specimens containing anticoagulants other than EDTA have not been validated for use with the **Determine HIV Early Detect** and may give incorrect results.
- Infants born to HIV-infected mothers may carry maternal antibodies and will test antibody positive until eighteen months of age, which may not necessarily indicate the true infection status of the new born. The use of HIV-1 p24 antigen testing to exclude infection in neonates (up to around eighteen months) is not recommended by CDC, because of poor sensitivity, especially in the presence of HIV antibody. Definitive diagnosis of HIV infection in early infancy requires other assays, including HIV nucleic acid test or viral culture.<sup>9</sup>
- Biotin treatment higher than 20mg per day may lead to decreased Ag bar intensity. Biotin concentrations up to 200 ng/mL in serum or plasma did not impact the sensitivity. There is no impact to Ab bar by biotin.

## PERFORMANCE CHARACTERISTICS

### SENSITIVITY

Sensitivity was evaluated by testing confirmed HIV antibody positive specimens, commercial seroconversion panels and HIV-1 subtype virus panels.

### 1. HIV Antibody Positive Specimens:

**Table I**  
**HIV Antibody Positive Specimens**

Types	Number of Specimens Tested	Reactive by Determine HIV Early Detect	Sensitivity
HIV-1	422	422	100.0%
HIV-1 non B subtypes *	56	56	100.0%
HIV-1 group O	4	4	100.0%
HIV-2	100	100	100.0%
HIV**	100	100	100.0%
<b>Total</b>	<b>682</b>	<b>682</b>	<b>100.0%</b>

\* Subtypes: A, C, D, F, G, H, J, K, and CRF01-AE, CRF02-AG, CRF03-AB, CRF05-DF, CRF09-A/U, CRF11-cpx

\*\* Specimens collected in Europe

A total of 682 confirmed HIV positive specimens were tested (Table I). The antibody diagnostic sensitivity of **Determine HIV Early Detect** on this population of specimens was calculated to be 100%.

### 2. HIV Seroconversion Panels: (See Table II on the last page.)

The sensitivity of the **Determine HIV Early Detect** was evaluated using 31 sets of seroconversion panels; each including early seroconversion panel members (98 early samples in total). The results were compared with the results of commercially available CE marked HIV 4<sup>th</sup> generation rapid immunochromatographic test kit, CE marked HIV 3<sup>rd</sup> generation rapid immunochromatographic test kit, HIV antibody EIA kit, HIV antigen EIA kit. Additionally, the results of 20 seroconversion panels were compared with the data of CE marked HIV 4<sup>th</sup> generation kit (Chemiluminescent microparticle immunoassay (CMIA)). On 28 of 31 panels, **Determine HIV Early Detect** detected HIV infection earlier than a 3<sup>rd</sup> generation rapid immunochromatographic test kit. On 11 of 31 panels, **Determine HIV Early Detect** detected HIV infection earlier than a rapid 4<sup>th</sup> generation immunochromatographic test kit.

Detection of six of 20 panels by **Determine HIV Early Detect** was delayed by 1 bleed date when compared to a 4<sup>th</sup> generation kit (CMIA).

Detection of eight of 31 panels by **Determine HIV Early Detect** was delayed by 1 bleed date when compared to HIV Antigen Kit (EIA).

### 3. HIV-1 p24 Subtype Antigens:

50 HIV virus panels were prepared by diluting 50 different cultured HIV isolates including HIV-1 group M subtypes: A(4), B(8), C(7), D(5), F(5), G(2), H(1), and CRF01-AE(10), CRF02-AG(2), HIV-1 group O(4) and HIV-2(2) with HIV negative human serum. The **Determine HIV Early Detect** detected all HIV-1 panels and did not detect HIV-2 panels.

### 4. Analytical Sensitivity of HIV-1 p24 antigen

The analytical sensitivity of the **Determine HIV Early Detect** was evaluated by testing WHO International Standard HIV-1 p24 Antigen (NIBSC code 90/636). The results demonstrated that the test could detect a concentration of 2 IU/mL HIV-1 p24 antigen.

## SPECIFICITY

A total of 2469 confirmed negative serum or plasma specimens were tested with the **Determine HIV Early Detect** and the specificity was Determined for the antibody (Ab) bar and for the antigen (Ag) bar (Table III). The specificity was 99.96% for the antibody (Ab) bar and 99.76% for the antigen (Ag) bar.

**Table III**  
**Specificity of Determine HIV Early Detect**

Population	Number of Tested Specimens	Determine HIV Early Detect					
		Ab bar		Ag bar		Combined Ab and Ag bar	
		Reactive	Non-reactive	Reactive	Non-reactive	Reactive	Non-reactive
Seronegative specimens*	1749	0	1749	2	1747	2	1747
Hospitalized specimens	218	1	217	2	216	3	215
Pregnant women	206	0	206	0	206	0	206
Potentially cross-reacting specimens**	296	0	296	2	294	2	294
Total	2469	1	2468	6	2463	7	2462
Specificity (%)		99.96%		99.76%		99.72%	

\* Including specimens collected in Europe (300), USA (1299) and Africa (150)

\*\* Disease states other than HIV, and potentially interfering substances: Rheumatoid factor, antinuclear antibody, systemic lupus erythematosus, high cholesterol, high total protein, high IgM, human anti mouse IgG, other infections (HBV, HCV, HTLV I or II, CMV, Toxo IgG, Syphilis, HSV 1/2 and EBV), Flu vaccinated patients, and dialysis patient specimens

## SAMPLE TYPE

### Table IV

#### Antibody sensitivity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

No. of matched specimens tested	Type of Specimens and No. of Reactive by Determine HIV Early Detect (Ab detection)				Correlation between matrices
	Serum	Plasma	Whole Blood venipuncture	Whole Blood fingerstick*	
25	25	25	25	NT	100.0%
25	25	25	25	25	100.0%

NT: not tested \*EDTA capillary tubes

**Multiple (matched) specimens:** Seropositive specimens from a total of 50 individuals from Africa were tested with the **Determine HIV Early Detect** (Table IV). Multiple (matched) specimens, 50 serum specimens, 50 plasma specimens, 50 whole blood (venipuncture) specimens and 25 whole blood (fingerstick) specimens were obtained from the donors in Africa. The results obtained from all specimen matrices showed 100% correlation, demonstrating that **Determine HIV Early Detect** gives identical results for these types of specimen matrices.

### Table V

#### Specificity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

No. of matched specimens tested	Type of Specimens and No. of Non-reactive by Determine HIV Early Detect				Correlation between matrices
	Serum	Plasma	Whole Blood venipuncture	Whole Blood fingerstick*	
25	25	25	25	NT	100.0%
25	25	25	25	25	100.0%

NT: not tested \*EDTA capillary tubes

**Multiple (matched) specimens:** Seronegative specimens from a total of 50 individuals from Africa were tested with the **Determine HIV Early Detect** (Table V). Multiple (matched) specimens, 50 serum specimens, 50 plasma specimens, 50 whole blood (venipuncture) specimens and 25 whole blood (fingerstick) specimens were obtained from the donors in Africa. The results obtained from all specimen matrices showed 100% correlation, demonstrating that **Determine HIV Early Detect** gives identical results for these types of specimen matrices.

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\* Vendor certificate data, \*\*PRB949 No.5 was not tested. NT: not tested

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**BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/  
BIBLIOGRAFIA**

1. Pilcher C, Eron JJ, Galvin S, et al. (2004) Acute HIV revisited: new opportunities for treatment and prevention. *The Journal of Clinical Investigations* 113(7): 937-945.
2. Respass RA, Rayfield MA and Dondero TJ (2001) Laboratory testing and rapid HIV assays: applications for HIV surveillance in hard-to-reach populations. *AIDS* 15 Supplement 3: S49-S59.
3. Louie B, Wong E, Klausner JD, et al. (2008) Assessment of rapid tests for detection of Human Immunodeficiency Virus-specific antibodies in recently infected individuals. *Journal of Clinical Microbiology* 46(4): 1494-1497.
4. Lyons MS, Lindsell CJ, Hawkins DA, et al. (2008) Contributions to early HIV diagnosis among patients linked to care vary by testing venue. *BMC Public Health* 8:220
5. Gürler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3):1581-1585.
6. Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
7. Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
8. EPA Guide for Infection Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
9. CDC: 2008 Case Definition; Human Immunodeficiency Virus Infection.
10. Delaney KP, Branson BM, Uniyal A, et al. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
11. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
12. O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.