Bioline™ HBsAg

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ENGLISH

About the test
[Intended use] The Bioline™ HBsAg test is an *in vitro* immunochromatographic assay designed for qualitative determination of HBsAg in human serum or plasma. Bioline™ HBsAg is intended only for professional use as an aid to diagnosis. Reactive specimens should be confirmed by a supplemental assay. This kit is for *in vitro* diagnostic use only. This product suitable for use on pregnant women, but in case of neonates, we haven't yet to perform the test. This test may not be suitable for diagnosis of early infection or blood donation screening.
[Test principle] This test device contains a membrane strip, which is pre-coated with mouse monoclonal anti-HBs antibody on test band region. The mouse monoclonal anti-HBs-colloid gold conjugate and the specimen move along the membrane chromatographically to the test region (T) and forms a visible line as the antibody-antigen-antibody gold particle complex forms.

- Materials provided**
The Bioline™ HBsAg test kit contains the following items to perform the assay:
- Bioline™ HBsAg [30 Tests/kit] Cat. No. 01FK10 :
 - 30 Test devices with desiccant in individual foil pouches
 - 1 Instructions for use
 - Bioline™ HBsAg [100 Tests/kit] Cat. No. 01FK11 :
 - 10 Multi-test devices with desiccant in individual foil pouch x 10
 - 1 Instructions for use


- Materials required but not provided**
- Micropipette

- Kit storage and stability**
1. The Bioline™ HBsAg should be stored between 1°C and 30°C. Do not freeze the kit.
 2. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
 3. Do not use the test kit beyond its expiration date.
 4. The shelf-life of the kits is as indicated on the outer package.
 5. Do not use the test kit if the pouch is damaged or the seal is broken.

- Specimen collection and storage**
1. [Plasma] Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture and then centrifuge blood to get plasma specimen. Anticoagulants such as heparin, EDTA, and citrate do not affect the test result.
 2. [Serum] Collect the whole blood into the collection tube (NOT containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture, leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
 3. If plasma or serum specimens are not tested immediately, they should be refrigerated at 2~8°C. For storage period longer than 3 days, freezing(below -20°C) is recommended. They should be brought to 15~30°C prior to use.
 4. Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

- Warnings**
1. For *in vitro* diagnostic use only. Do not reuse the test device.
 2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
 3. Do not eat or smoke while handling specimens.
 4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
 5. Avoid splashing or aerosol formation of specimen.
 6. Clean up spills thoroughly using an appropriate disinfectant.
 7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.

- Test procedure**
1. Allow all kit components and specimen to room temperature (between 15°C and 30°C) prior to testing.
 2. Remove the test device from the foil pouch and place it on a flat, dry surface.
 3. Dispense 100µl of plasma or serum specimen into the specimen well (Figure 1).
 4. As the test begins to work, you will see purple color move across the result window in the center of the test device.
 5. Interpret test results at 20 minutes.

◆ SP 100 µl : Serum 100µl or Plasma 100µl

- Test interpretation**
1. A purple control band will appear in the left section of the results window to show that the test is working properly.
 2. The right section of the results window indicates the test results. If another purple band appears in the right section of the results window, this is the test band.
- **Negative Result** : The presence of only control band (C) within the results window indicates a negative result.
-
- **Positive Result** : The presence of the test line (T) and the control line (C) within the results window, regardless of which band appears first, indicates a positive result.
-
- **Caution:** The presence of any test line, no matter how faint, the result is considered positive.
-
-
-

Test limitations
A negative result does not preclude the possibility of infection with HBV. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Internal quality control
The Bioline™ HBsAg test has a test line (T) and a control line (C) on the surface of the device. Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the active ingredients of main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Origin	HBsAg positive specimens	HBsAg negative specimens	Total
Africa	16	43	59
Asia	30	30	60
Europe	30	67	97
Latin America	22	38	60
Total	98	178	276

In this multi-site evaluation of 276 specimens, we found the relative sensitivity is 100% (98/98), the relative specificity is 100% (178/178). The results are summarized in the following tables.

Results		Bioline™ HBsAg			95% CI
		Positive	Negative	Total Results	
Reference	Positive (98)	98	0	98	96.3 - 100%
	Negative (178)	0	178	178	97.9 - 100%
Total		98	178	276	

Product Disclaimer:
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning:
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect positive or negative diagnosis using this product.

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© 2021 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.
Abbott Diagnostics Korea Inc.