Abbott

Bioline®

HIV 1/2 3.0

The 3rd Generation of antibodies to HIV-1/HIV-2 Test Le 3^e de g n ration d'anticorps anti-VIH-1/VIH-2 The 3^{er} g n ration d'anticorps pour la prueba de VIH 1 o 2 A 3^a g n ration de anticorpos do teste VIH-1/VIH-2

ENGLISH

About the test

[Introduction] Le VIH (Human Immunodeficiency Virus) is recognized as the etiologic agent of Acquired Immune Deficiency Syndrome (AIDS). The virus is transmitted by sexual contact, exposure to infected blood, certain body fluids or tissues, and from mother to fetus or child during the perinatal period. HIV-1 has been isolated from patients with AIDS and AIDS related complex, and from healthy persons with high potential risk of developing AIDS. Patients with HIV-2 are found primarily in parts of West Africa. Its course is marked by increasing levels of viral replication and the emergence of more virulent viral strains. HIV-1 and HIV-2 are similar in their morphology, cell targets, host interaction and genetic structure. Serological studies have determined that HIV-1 and HIV-2 have multiple common epitopes in core antigens but much less in the envelope antigens. The clinical diagnosis of HIV may include the detection of antibodies to HIV-2 in human plasma or serum by immunosorbent assay. The presence of HIV can be identified by detection of antibodies to HIV-1 in human serum, plasma and whole blood by immunoassay. This advanced assay utilizes recombinant antigen targeted against immunogenic proteins. The major immunogenic antigens of these proteins are HIV-gp120 and HIV-gp41 and gp36.

[Test principle] The Bioline® HIV 1/2 3.0 test contains a membrane strip, which is pre-coated with recombinant HIV-1 capture antigen (gp41, p24) on test line 1 region and with recombinant HIV-2 capture antigen (gp41, p24) on test line 2 region respectively. The recombinant HIV-1/2 antigen (gp41, p24 and gp36)-coiled gold conjugate and the sample move along the membrane chromatographically to the test region (T) and forms a visible line as the antigen-antibody-antigen complex. The test results are visible within 10 minutes.

[Intended use] The Bioline® HIV 1/2 3.0 is a rapid, qualitative test for the detection of antibodies to all isotopes (IgG, IgM, IgA) specific to HIV-1 and HIV-2 simultaneously in human serum, plasma or whole blood. The Bioline® HIV 1/2 3.0 kit is intended only for professional use and for *in vitro* diagnostic use. This test may not be suitable for diagnosis of early infection or blood donation screening. Positive samples should be confirmed by a supplemental assay such as ELISA or Western Blot Test. The performance of the assay has not been established for populations of infants, children or adults.

Materials provided and active ingredients of main components

- The Bioline® HIV 1/2 3.0 test kit contains the following items to perform the assay:
 - 1-1. Bioline® HIV 1/2 3.0 [30 Tests/Kit] Cat. No. 03FK10;
 - 30 Test devices with desiccant in individual foil pouches.
 - Assay diluent (1 x 4 mL/vial).
 - 1 Instructions for use.
- 2-1. Bioline® HIV 1/2 3.0 [25 Tests/Kit] Cat. No. 03FK16;
- 25 Test devices with desiccant in individual foil pouches.
- Assay diluent (1 x 4 mL/vial).
- 25 Capillary pipettes (20 µl), 25 Sterile lancets, 25 Alcohol swabs.
- 1 Instructions for use.
- Active ingredients of main components
 - 1-1. Test strip including Gold conjugate: Recombinant HIV-1 gp41, p24 antigen, HIV-2 gp36 antigen, colloidal gold (10.0±0.2 µg), Test line 1: Recombinant HIV-1 antigen (gp41, p24) (0.025±0.025 µg), Test line 2: Recombinant HIV-2 antigen (gp36) (0.05±0.1 µg), Control line: Anti-HIV serum (0.75±0.05 µg).
 - Assay diluent: 50mM Tris-HCl Buffer, Sodium azide (s.a.s.)

Materials required but not provided

- Micropipette, Protective gloves, Timer, Biohazard container

Kit storage and stability

- The test kit should be stored at a temperature between 1°C and 30°C. Do not freeze the kit or its components.
- Assay diluent may be opened and respsted for each assay. Cap should be firmly sealed between each use. Assay diluent is stable until expiration date if kept at 1 - 30°C.
- The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from foil pouch.
- Do not use the test kit beyond its expiration date. The shelf life of the kit is indicated on the outer package.
- Do not use test kit if the pouch is damaged or the seal is broken.

Warnings

- The test devices are for *in vitro* diagnostic use only. Do not reuse the test device.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be proficient.
- Do not use the pipette by mouth, smoke, drink, eat, apply cosmetics, or handle contact lenses in areas where gloves are being handled.
- Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.
- Do not mix or interchange different specimens.
- Do not use the test device for anything other than the intended purpose.
- Avoid splashover or aerosol formation of specimen and assay diluent.
- Do not mix or interchange components among different lots or those for other products.
- Do not drink assay diluent.
- Care should be taken to avoid contamination of the bottle nozzle when dropping assay diluent into the specimen well.
- The assay diluent contains a proprietary antimicrobial agent, sodium azide, which prevents any growth of bacteria. Do not use the assay diluent for anything other than the intended purpose. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.
- The assay diluent contains sodium azide, which may react with lead or copper plating to form highly explosive metal azide compounds. When disposing of these reagents through plumbing/filters, flush with a large volume of water to prevent azide build-up in drains. Safety data sheet available for professional user on request.

Specimen collection and handling

- Whole blood**
 - [Collection by venipuncture]** Using vacuumtainer, collect whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
 - If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8°C.
 - If stored at 2 - 8°C, the blood specimen must be tested within 3 days of refrigeration.
 - Do not use a blood specimen to room for more than 3 days; it can cause a nonspecific reaction.
 - Bring blood specimens to room temperature (15 - 30°C) prior to use.
- Plasma or serum**
 - [Plasma]** Using vacuumturer, draw the whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) and then centrifuge blood to obtain the plasma specimen.
 - [Serum]** Using vacuumturer, draw the whole blood into the collection tube (NOT containing anticoagulants) then leave for 30 minutes to allow blood coagulation to occur. Centrifuge the tube to separate the serum and discard the clot.
 - If plasma or serum specimens are not tested immediately, they must be refrigerated at 2 - 8°C. For storage period longer than 2 weeks, freezing (below -20°C) is required. Bring plasma or serum specimens to room temperature (15 - 30°C) prior to use.
 - Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified by centrifugation prior to assaying.
- Precautions**
 - Anticoagulants including heparin, EDTA and sodium citrate do not affect the test result.
 - Use of other anticoagulants has not been validated. Their use may affect the test result.
 - Use new pipette tips for each specimen in order to avoid cross-contamination of other specimens which could cause erroneous results.

Glossary of symbols / Glossaire des symboles / Glosario de s mbolos / Gloss rio de s mbolos

	Temperature limitation Limitation de temperature L�mite de temperatura Limite de temperatura		Lot Number Limitation de temperature N�mero de lote Numero de lote		Use by / Exp Date de pr�paration / Utiliser avant le Fecha de vencimiento / Fecha Prazo de validade / Exp
	In vitro diagnostic medical device Dispositif m�dical de diagnostic Dispositivo m�dico de diagn�stico in vitro Dispositivo m�dico de diagn�stico in vitro		Catalog Number N�mero de catalogue Numero de catalogo		Consult instructions for use Consulte las instrucciones Consulte las instrucciones de uso Consulte as instru�es de utiliza��o
	Do not reuse No reutilizar N� reutilizar		Contains sufficient for one tests Contenu suffisant pour uno tests Cont�nido suficiente para realizar uno prubas Cont�nido suficiente para uno testes		Keep away from sunlight Conservar a l�o da luz da luzimre do sol Mantener alejado de la luz solar Manten afastado da luz solar
	Keep dry Garder au sec Mantener seco Mantenr seco		Biological risks Riscos biol�gicos Riscos biol�gicos Riscos biol�gicos		Do not use if package is damaged Ne pas utiliser si l'emballage est endommag�e No utilizar si el paquete est da�ado N� utilizar se o embalagem est�v danificada
	Manufacturer Fabricant Fabricante Fabricante		Date of manufacture Date de fabrication Data de fabrica��o Data de fabrica��o		Caution Attention Precau��o Aten��o

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PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

- 1** **EN** Now, open the package and look for the following:

 1. Test device with desiccant in individual foil pouch
 2. Assay diluent
 3. Instructions for use

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Incluído apenas para o n.º de catálogo 03FK16

 4. Pipeta capilar (20 µl)
 5. Lanceta estéril
 6. Zaragatoa com álcool

- 2** **EN** Carefully read the instructions on how to use the Bioline™ HIV 1/2 3.0 test kit.

ES Lea atentamente las instrucciones sobre cómo utilizar el kit de prueba Bioline™ HIV 1/2 3.0.
- FR** Lisez attentivement les instructions d'utilisation du kit de test Bioline™ HIV 1/2 3.0.

PT Leia atentamente as instruções sobre como utilizar o kit de teste Bioline™ HIV 1/2 3.0.
- 3** **EN** Look at the expiration date at the back of the foil pouch. If the expiration date has passed, use another kit. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

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- 4** **EN** Open the foil pouch and look for the following:

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
- 1** **EN** Dispense 10 µl of plasma or serum specimen or 20 µl of whole blood specimen into the specimen well marked "S"

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- 2** **EN** Dispense 4 drops (approximately 120 µl) of assay diluent into the specimen well "S". Hold bottle vertically while dispensing. Do not let bottle nozzle touch device in order to avoid cross-contamination.

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- 3** **EN** Interpret the result 10 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 10 min or after 20 min) may provide in false results.

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- 1** **EN** Clean the area to be lanced with an alcohol swab.

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- 3** **EN** Immerse the open end of a new capillary pipette (20 µl) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

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- 1** **EN** Dispense 20 µl of drawn whole blood specimen in the specimen well marked "S".

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INTERPRETATION / INTERPRÉTATION / INTERPRETACIÓN / INTERPRETAÇÃO

Negative / Négatif / Negativo

- EN** The presence of only the control line (C) within the result window indicates a negative result.

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Positive / Positif / Positivo

- EN** ⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered positive.

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HIV-1 Positive / Positif / Positivo

- 2 LINES / 2 LIGNES / 2 LÍNEAS / 2 LINHAS**

EN The presence of both test line 1 (T) and the control line (C) indicates a positive result for anti-HIV-1.

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Weak / Faibles / Tenue / Ténue
-

3 LINES / 3 LIGNES / 3 LÍNEAS / 3 LINHAS

- EN** If the color intensity of the test line 1 is darker than one of test line 2 in the result window, you can interpret the result as anti-HIV-1 positive.

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HIV-2 Positive / Positif / Positivo

- 2 LINES / 2 LIGNES / 2 LÍNEAS / 2 LINHAS**

EN The presence of both test line 2 (T) and the control line (C) indicates a positive result for anti-HIV-2.

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Weak / Faibles / Tenue / Ténue
-

3 LINES / 3 LIGNES / 3 LÍNEAS / 3 LINHAS

- EN** If the color intensity of the test line 2 is darker than one of test line 1 in the result window, you can interpret the result as anti-HIV-2 positive.

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- EN** [Both HIV-1 and HIV-2 Positive]
The presence of the test line 1 (T) and the test line 2 (T) and the control line (C) indicates a positive result for anti-HIV-1 and/or anti-HIV-2.

⚠ **Caution:** Although a positive result for anti-HIV-1 and anti-HIV-2 in one patient is a rare case, it's possible as there is an homology in the amino acid sequence between anti-HIV-1 and anti-HIV-2. To determine the virus type or diagnose a co-infection accurately, you must perform a confirmatory test such as Western Blot etc.
- FR** [Positifs aux deux VIH-1 et VIH-2]
La présence de la ligne de test 1 (T), de la ligne de test 2 (T) et de la ligne de contrôle (C) indique un résultat positif pour les anticorps anti-VIH-1 et/ou anti-VIH-2.

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Invalid / Non valide / No válido / Inválido

- EN** Absence of the control line (C) and/or presence of a pink/purple smear in the result window indicates an invalid result. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

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