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1. D. T. Selewski and J. M. Symons, "Acute kidney injury," *Pediatr. Rev.*, vol. 35, no. 1, pp. 30–41, 2014.
2. L. Koppe, D. Fouque, and C. O. Soulage, "Metabolic Abnormalities in Diabetes and Kidney Disease: Role of Uremic Toxins," *Curr. Diab. Rep.*, vol. 18, no. 10, 2018.
3. K. Nakamura, J. Kido, H. Mitsubuchi, and F. Endo, "Diagnosis and treatment of urea cycle disorder in Japan.," *Pediatr. Int.*, vol. 56, no. 4, pp. 506–509, 2014.
4. M. Ostermann, "Diagnosis of acute kidney injury: Kidney Disease Improving Global Outcomes criteria and beyond.," *Curr. Opin. Crit. Care*, vol. 20, no. 6, pp. 581–587, 2014.
5. A. H. B. Wu, *Tietz clinical guide to Laboratory tests*, 4th Ed. Saunders Elsevier, 2006.
6. Kaplan A. Urea. Kaplan A et al. *Clin Chem The C.V. Mosby Co.* St Louis. Toronto. Princeton 1984; 1257-1260 and 437 and 418.
7. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
8. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.

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Ref: 1001331

Cont.

R1: 2 x 150 mL; R2: 2 x 150 mL; R3: 2 → 150 mL
CAL: 1 x 5 mL

Ref: 1001329

R1: 5 x 50 mL; R2: 5 x 50 mL; R3: 5 → 50 mL,
CAL: 1 x 5 mL

Quantitative determination of urea



IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD

Urea in the sample is hydrolyzed enzymatically into ammonia (NH₄⁺) and carbon dioxide (CO₂).

Ammonia ions formed reacts with salicylate and hypochlorite (NaClO), in presence of the catalyst nitroprusside, to form a green indophenol:

The intensity of the color formed is proportional to the urea concentration in the sample⁶.

CLINICAL SIGNIFICANCE

Urea is the final result of the metabolism of proteins; it is formed in the liver from its destruction.

Elevated urea can appear in blood (uremia) in: diets with excess of proteins, renal diseases, heart failure, gastrointestinal hemorrhage, dehydration or renal obstruction⁶.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1 Buffer	Phosphate pH 6,7	50 mmol/L
	EDTA	2 mmol/L
	Sodium salicylate	400 mmol/L
	Sodium nitroprusside	10 mmol/L
R 2 NaClO	Sodium hypochlorite (NaClO)	140 mmol/L
	Sodium hydroxide (NaOH)	150 mmol/L
R 3 Enzymes	Urease	62500 U/L
	UREA CAL	Urea aqueous primary standard
		50 mg/dL

PRECAUTIONS

R2: H314-Causes severe skin burns and eye damage. H412-Harmful to aquatic life with long lasting effects.

Follow the precautionary statements given in MSDS and label of the product.

PREPARATION

- Working reagent (WR): Dissolve (→) one tablet R 3 Enzymes in one bottle of R 1 Buffer. Cap and mix gently to dissolve contents.

Stability: 4 weeks in the refrigerator (2-8°C) or 7 days at room temperature (15-25°C).

- R 2 NaClO is ready to use.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use.

Do not use reagents over the expiration date.

Do not use the tablets if appears broken.

Signs of reagent deterioration:

- Presence of particles and turbidity.
- Blank absorbance (A) at 580 nm ≥ 0,32.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 580 nm.
- Matched cuvettes 1,0 cm light path.
- General laboratory equipment (Note 2).

SAMPLES

- Serum or heparinized / EDTA plasma^{1,4}.

- Urine³: Dilute sample 1/50 in distilled water. Mix.

In all described sample types, urea is stable for at least 7 days at 2 - 8 °C or 35 days at -20 °C.

PROCEDURE

1. Assay conditions:

Wavelength: 580 nm

Cuvette: 1 cm light path

Temperature: 37°C / 15-25°C

2. Adjust the instrument to zero with distilled water.

3. Pipette into a cuvette: (Note 4)

	Blank	Standard	Sample
WR (mL)	1,0	1,0	1,0
Standard ^(Note 1,3) (μL)	--	10	--
Sample (μL)	--	--	10

4. Mix and incubate 5 min at 37°C or 10 min at room temperature (15-25°C).

5. Pipette:

	Blank	Standard	Sample
R 2 (mL)	1,0	1,0	1,0

6. Mix and incubate 5 min at 37°C or 10 min at room temperature (15-25°C).

7. Read the absorbance (A) of the samples and calibrator, against the Blank. The colour is stable for at least 30 minutes at 15-25°C.

CALCULATIONS

(A) Sample – (A) Blank x 50 (Standard conc.) = mg/dL urea in the sample

(A) Standard – (A) Blank

mg/dL Urea x 0,466 = mg/dL Urea BUN (Blood Urea Nitrogen)⁶.

Conversion factor: mg/dL x 0,1665 = mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINTROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES⁵

Serum : 12.6 - 42.6 mg/dL (2.1 - 7.1 mmol/L)

Urine : 25.8 - 42.6 gr/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

1. **Trueness:** Trueness was assessed by a method comparison with real samples. The obtained results can be consulted in the upcoming section titled.

2. **Precision:** repeatability and reproducibility were evaluated based on the CLSI protocol EP05-A3:

N	Mean (mg/dL)	Repeatability	
		SD	CV (%)
80	40.4	0.41	1.0
	113	0.23	0.2

N	Mean (mg/dL)	Reproducibility	
		SD	CV (%)
90	43.8	2.64	6.0
	117	5.84	5.0

3. **Analytical sensitivity:** the observed analytical sensitivity was 0.00588 Abs/mg/dL at 580 nm.

4. **Analytical specificity:** bilirubin (21.2 mg/dL), triglycerides (4097 mg/dL), hemoglobin (312 mg/dL) do not interfere. Do not use ammonium salts as anticoagulant. A list of drugs and other interfering substances with urea determination has been reported by Young et. al^{7,8}.

5. **Limit of detection (LoD):** values less than 0.27 mg/dL give non-reproducible results.

6. **Measurement range:** from quantitation limit of 1.50 mg/dL to linearity limit of 250 mg/dL. If the results obtained are greater than the linearity limit, then dilute the sample 1/2 with NaCl 9 g/L and retested again.

7. **Method comparison:** results obtained using this reagent (y) were compared to those obtained using a commercial reagent (x) with similar characteristics. 50 samples of different urea concentrations were assayed. The correlation coefficient (r) and the regression equation were:

y = 1.081 + 0.973x r = 0.999

The results of the performance characteristics depend on the analyser used.

NOTES

1. UREA CAL: Proceed carefully with this product because due its nature it can get contaminated easily.
2. Glassware and distilled water must be free of ammonia and ammonium salts⁶.
3. Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.
4. Use clean disposable pipette tips for its dispensation.
5. **SPINREACT has instruction sheets for several automatic analyzers.**

BIBLIOGRAPHY

1. D. T. Selewski and J. M. Symons, "Acute kidney injury," *Pediatr. Rev.*, vol. 35, no. 1, pp. 30–41, 2014.
2. L. Koppe, D. Fouque, and C. O. Soulage, "Metabolic Abnormalities in Diabetes and Kidney Disease: Role of Uremic Toxins," *Curr. Diab. Rep.*, vol. 18, no. 10, 2018.
3. K. Nakamura, J. Kido, H. Mitsubuchi, and F. Endo, "Diagnosis and treatment of urea cycle disorder in Japan," *Pediatr. Int.*, vol. 56, no. 4, pp. 506–509, 2014.
4. M. Ostermann, "Diagnosis of acute kidney injury: Kidney Disease Improving Global Outcomes criteria and beyond," *Curr. Opin. Crit. Care*, vol. 20, no. 6, pp. 581–587, 2014.
5. A. H. B. Wu, *Tietz clinical guide to Laboratory tests*, 4th Ed. Saunders Elsevier, 2006.
6. Kaplan A. Urea. Kaplan A et al. *Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton* 1984; 1257-1260 and 437 and 418.
7. Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
8. Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.

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BIBLIOGRAPHIE

- D. T. Selewski and J. M. Symons, "Acute kidney injury," *Pediatr. Rev.*, vol. 35, no. 1, pp. 30-41, 2014.
- L. Koppe, D. Fouque, and C. O. Soulage, "Metabolic Abnormalities in Diabetes and Kidney Disease: Role of Uremic Toxins," *Curr. Diab. Rep.*, vol. 18, no. 10, 2018.
- K. Nakamura, J. Kido, H. Mitsubuchi, and F. Endo, "Diagnosis and treatment of urea cycle disorder in Japan," *Pediatr. Int.*, vol. 56, no. 4, pp. 506-509, 2014.
- M. Ostermann, "Diagnosis of acute kidney injury: Kidney Disease Improving Global Outcomes criteria and beyond," *Curr. Opin. Crit. Care*, vol. 20, no. 6, pp. 581-587, 2014.
- A. H. B. Wu, *Tietz clinical guide to Laboratory tests*, 4th Ed. Saunders Elsevier, 2006.
- Kaplan A. Urea. Kaplan A et al. *Clin Chem The C.V. Mosby Co.* St Louis. Toronto. Princeton 1984; 1257-1260 and 437 and 418.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.

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