# LDL Colesterol

## Reactivo Precipitante

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# LDL Cholesterol

## *Reactivo Precipitante*

For the separation of the low density lipoproteins  
(LDL) in serum

### SUMMARY

The approximate cholesterol content in each lipoprotein family is (in % per weight unit): 1% in chylomicrons, 18% in VLDL, 50% in LDL and 23% in HDL. Each family has a different biological activity, thus the clinical significance of a cholesterol increase depends on the lipoprotein or lipoproteins that are in excess.

The regulatory mechanisms of the plasma lipoproteins levels are very complex and can be affected by several factors (genetic, environmental, physiological or pathological). It is possible to find total cholesterol values near the normal range along with alterations in the lipoprotein fractions.

The HDL and LDL have been the most studied due to their important biological activity:

- the LDL, product of the VLDL metabolism in plasma, are the ones that transport the exogenous cholesterol (and in much less proportion, endogenous) to the interior of the cells;
- the HDL, synthesized in the liver, remove the cholesterol that is not used by the cells (within certain concentration limits), transporting it to the liver for its degradation.

Several epidemiological studies have confirmed that the LDL cholesterol excess in reference to a critical value (1.9 g/l) must be considered as a risk factor for the development of coronary heart disease, considering that the HDL's protective effect only seems to have relevance within certain concentrations range of circulating cholesterol. These findings allow to deduce that isolated HDL or LDL values cannot be taken as predictive risk indexes, but it is necessary to make a lipid profile with the total cholesterol, HDL cholesterol and LDL cholesterol values.

### PRINCIPLE

Low-density lipoproteins (LDL or  $\beta$ -lipoproteins) are separated from the serum selectively precipitating them by adding high molecular weight polymers. In the supernatant, after centrifugation, the other lipoproteins (HDL and VLDL) remain; the determination of cholesterol associated to them is performed using the enzymatic system cholesterol-oxidase/peroxidase with Trinder colorimetry (Phenol/4-AP). The difference between the total cholesterol and the one determined in the supernatant will give the cholesterol associated to the LDL.

### PROVIDED REAGENTS

**A. Reagent A** (Precipitating Reagent): 1 g/l solution of polyvinyl sulfate dissolved in polyetyleneglycol (MW: 600) at 25%, pH 6.7.

### NON-PROVIDED REAGENTS

Wiener lab.'s **Colestat enzimático**, **Colestat enzimático AA** or **Colestat enzimático AA líquida**.

### INSTRUCTIONS FOR USE

**Reagent A:** ready to use.

### WARNINGS

The Reagent is for "in vitro" diagnostic use.

Use the reagents according to the working procedures for clinical laboratories.

The reagents and samples should be discarded according to the local regulations in force.

### STABILITY AND STORAGE INSTRUCTIONS

**Reagent A:** stable in refrigerator (2-10°C) until the expiration date shown on the box.

### INSTABILITY OR DETERIORATION OF REAGENTS

Any change in the coloration or in other physical feature of the reagent may be a sign of its deterioration.

### SAMPLE

Serum

**a) Collection:** the patient must be on a 12 to 16-hour fasting. Obtain the serum in the usual way and separate the clot within an hour from the extraction.

**b) Additives:** not required.

**c) Known interfering substances:** hypertriglyceridemic sera (with chylomicronemia) produce turbid supernatants; bilirubin in levels above 50 mg/l, interferes.

**d) Stability and storage instructions:** serum should be preferably fresh. In case the test cannot be performed immediately, the sample can be kept refrigerated (2-10°C) for no longer than 24 hours from the extraction. Do not freeze.

### REQUIRED MATERIAL (non-provided)

- Centrifuge.
- Spectrophotometer or photocolormeter.
- Micropipettes and pipettes for measuring the stated volumes
- Kahn tubes.
- Spectrophotometric cuvettes.
- Water bath at 37°C.
- Stopwatch.

### PROCEDURE

In a Kahn tube place:

|           |             |
|-----------|-------------|
| Sample    | 200 $\mu$ l |
| Reagent A | 100 $\mu$ l |

Homogenize by shaking (not by inversion) for 20 seconds and place in a water bath at 20-25°C for 15 minutes. Centrifuge for 15 minutes at 3000 rpm. Immediately separate the supernatant. See PROCEDURE LIMITATIONS.

Use the supernatant as Sample for the colorimetric assay. In three photocolormeter tubes labeled B (Blank), S (Standard) and U (Unknown) place:

|                        | B    | S     | U      |
|------------------------|------|-------|--------|
| <b>Supernatant</b>     | -    | -     | 100 ul |
| <b>Standard</b>        | -    | 20 ul | -      |
| <b>Working Reagent</b> | 2 ml | 2 ml  | 2 ml   |

Mix and incubate for 5 minutes at 37°C, if the Working Reagent of **Colestat enzimático AA Línea Líquida** is used or 15 minutes at 37°C if **Colestat enzimático** Working Reagent is employed. Remove from the bath and let cool. Read at 505 nm in spectrophotometer or at 490-530 nm in photocolormeter with green filter, setting the instrument to zero with the Blank.

## CALCULATIONS

LDL Cholesterol (g/l) = Total Cholesterol (\*) - (U x f)

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$FV_E$  = extract final volume = 0.3 ml

$V_S$  = processed sample volume = 0.2 ml (200 ul)

$RV_E$  = reaction volume with extract = 2.1 ml

$RV_{Std}$  = reaction volume with Standard = 2.02 ml

$V_{Std}$  = Standard volume in the reaction = 0.02 ml (20 ul)

$V_E$  = extract volume in the reaction = 0.1 ml (100 ul)

If different volumes are used the 0.624 factor varies and must be calculated again.

## REFERENCE VALUES

The expert panel from the National Cholesterol Education Program (NCEP) provides the following LDL Cholesterol values related to the risk of acquiring coronary heart disease (CHD):

- **Low or no risk** (normal individuals): LDL Cholesterol values below 1.29 g/l.

- **Moderate or high risk** (individuals with probability of acquiring CHD): values between 1.30 and 1.89 mg/l.

- **Very High risk** (individuals suspected from suffering from CHD): LDL Cholesterol values  $\geq 1.90$  g/l.

However, each laboratory should establish its own references values.

## PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

When icteric samples are processed, the sera should be diluted 1:2 or 1:3 with saline solution, using the usual procedure, taking into account the dilution factor for the calculations.

When the LDL cholesterol cannot be completely separated in a common centrifuge, due to high levels of triglycerides, the determination can be performed in this manner:

Follow the instructions under PROCEDURE up to the incubation in water bath at 20-25°C for 15 minutes, then place the reaction mixture in capillary tubes and centrifuge in a microhematocrit centrifuge at 10,000 rpm for 5 minutes. Cut the tube discarding the precipitate and use the clean supernatant for the test.

In all the cases in which the cholesterol values associated to HDL and VLDL ( $U \times f$ , according to the calculation) are above 1 g/l, the determination must be repeated with the same sera using 100 ul sample and 200 ul Reagent A. Continue with the described technique, taking 100 ul supernatant and using the 1.248 factor for the calculations.

Do not let the precipitate in contact with the supernatant, as a redissolution may take place producing high reading values with the consequent decrease in the LDL Cholesterol values.

## PERFORMANCE

**a) Reproducibility:** processing replicates of one sample on the same day, the following results were obtained:

| Level    | S.D.           | C.V.  |
|----------|----------------|-------|
| 1.14 g/l | $\pm 0.03$ g/l | 2.6 % |
| 2.03 g/l | $\pm 0.04$ g/l | 2.0 % |

**b) Linearity:** the reaction is linear up to 5 g/l.

**c) Detection limit:** depends on the spectrophotometer used. For a 0.001 O.D. reading, the minimum visible change of concentration will be of approximately 0.0063 g/l.

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Wiener Laboratorios S.A.I.C.