Uric Acid

fast

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REFERÊNCIAS

- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, Sixth Edition. Elsevier 2018.
- Tietz Clinical Guide to Laboratory Tests, 4th edition. St. Louis (MO): Saunders Elsevier 2006.
- CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. Evaluation of Linearity of Quantitative Measurement Procedures. 2nd ed. CLSI guideline EP06. Clinical and Laboratory Standards Institute; 2020.
- CLSI. Interference Testing in Clinical Chemistry. 3rd ed. - CLSI guideline EP07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
- CLSI. Measurement Procedure Comparasen and Bias Estimation Using Patient Samples. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.

Uric Acid

fast

For determination of uric acid in serum, plasma or urine

INTENDED USE


Enzymatic method for the quantitative determination of uric acid in serum, plasma or urine in automated analyzers. Intended for professional use in clinical analysis laboratories, medical clinics, and hospitals.

SUMMARY

In humans, uric acid (2,6,8-trihydroxypurine) is the main product of purine nucleoside catabolism. It is eliminated mainly by urinary excretion. The most frequent clinical use for uric acid determination in serum or plasma is the assessment of gout risk and its therapeutic monitoring. Other clinical conditions where serum urate measurements have potential utility include numerous renal and metabolic disorders such as chronic renal failure, leukemia, psoriasis, preeclampsia, alcoholism, and in patients receiving cytotoxic drugs. Uric acid measurements in urine may play a role in evaluating the cause of hyperuricemia and assessing the risk of kidney stone formation.

PRINCIPLE

Enzyme uricase reacts with uric acid present in the sample to produce allantoin, carbon dioxide and hydrogen peroxide. Hydrogen peroxide, in a reaction catalyzed by the POD enzyme, produces the oxidative condensation of 4-AP with TOOS (N-ethyl-N-(2-hydroxy-3-sulfopropyl)-3-methylaniline), resulting in a colored quinoneimine that is monitored photometrically at 575 nm.

UOD: uricase

POD: peroxidase

4-AP: 4-aminophenazone

TOOS: N-ethyl-N-(2-hydroxy-3-sulfopropyl)-3-methylaniline

PROVIDED REAGENTS

S. Standard*: uric acid solution 10 mg/dl.

A. Reactivo A: solution containing Good buffers pH 7.8, N-ethyl-N-(2-hydroxy-3-sulfopropyl)-3-methylaniline (TOOS), 4-aminophenazone (4-AP), uricase (UOD), peroxidase (POD), and potassium ferrocyanide.

Good Buffers.....50 mmol/l

UOD.....≥ 150 U/l

POD.....≥ 1000 U/l

4-AP.....0.10 mmol/l

Potassium ferrocyanide.....6 μmol/l
TOOS.....≤ 2.0 mmol/l

NON-PROVIDED REAGENTS

- Wiener lab **Calibrador A plus**.
- Distilled water.

INSTRUCTIONS FOR USE

Provided Reagents: ready to use.

WARNINGS

The reagents are for *in vitro* diagnostic use.

Do not pipette by mouth. Avoid contact with skin and eyes. In case of spillage or splashes, wash the affected area with plenty of water.

Use the reagents following the usual precautions for clinical chemistry laboratory work. All reagents and samples should be handled according to local regulations.

All reagents and samples should be discarded according to local regulations.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable in refrigerator (2-10°C) until the expiration date indicated on the box. Once opened, they should not remain uncovered and out of the refrigerator for long periods of time. Avoid contamination.

INSTABILITY OR DETERIORATION OF REAGENTS

Difficulty in obtaining control values within the assigned range (e.g. **Standatrol S-E 2 niveles**) is an indication of Reagent deterioration. In such case, discard. Turbidity is an indication of Reagent deterioration.

Discard when Target readings are > 0.200 O.D. or Standard readings are abnormally low.

SAMPLE

Serum, plasma or urine

a) Collection: Serum or plasma should be obtained in the usual manner. Separate the clot as soon as possible, within two hours after collection. If the sample is urine, preferably use fresh urine.

b) Additives: heparin or EDTA when plasma is used as sample.

c) Known interfering substances:

- Medications: Reducing substances such as ascorbic acid (vitamin C), Buscopan (Hyoscine butylbromide), etc. in high doses interfere. For this reason, medication should be discontinued, whenever possible, 24 hours prior to sampling.

- No interferences are observed for bilirubin up to 7.5 mg/dL, triglycerides up to 300 mg/dL, haemoglobin up to 740 mg/L. Refer to Young's bibliography for the effects of drugs in the present method.

d) Stability and storage instructions: centrifuge and separate the serum from the clot within two hours after extraction. If samples are not processed immediately, they can be stored for 1 week in refrigerator (2-10°C) and for 30 days at -20°C. Urine samples can be stored for 4 days at 20-25°C at pH > 8. Do not refrigerate or freeze.

REQUIRED MATERIAL (not provided)

- Spectrophotometer or photocolormeter.
- Suitable volumetric material.
- Parallel-sided spectrophotometric tubes or cuvettes.
- Water bath at 37°C.
- Clock or timer.

REACTION CONDITIONS

- Wavelength: 570 nm in spectrophotometer.
- Reaction temperature: 37°C
- Reaction time: 100 seconds
- Sample volume: 20 µL
- Final reaction volume: 900 µL

Sample and Reagent volumes can be decreased or increased proportionally (e.g. 50 µL Sample + 2250 µL Reagent A or 10 µL + 450 µL).

PROCEDURE

I- SERUM OR PLASMA TECHNIQUE

In three tubes or spectrophotometric cuvettes marked B (Blank), S (Standard or Calibrator) and U (Unknown), place:

	B	S	U
Demineralized water	20 µL	-	-
Standard or Calibrator	-	20 µL	-
Sample	.	-	20 µL
Reagent A	900 µL	900 µL	900 µL

Mix gently and incubate for 100 seconds in water bath at 37°C. Remove and read in spectrophotometer at 570 nm, taking the instrument to zero with the blank.

II- URINE TECHNIQUE

Use technique I diluting the urine 1/10 with water or saline solution. To calculate the results, multiply by the dilution factor used.

CALCULATION OF RESULTS

Uric acid (mg/L) = U x f

$$f = \frac{10 \text{ mg/dL}}{S}$$

U: Absorbance of Sample – absorbance of Blank
S: Absorbance of Standard – absorbance of Blank

QUALITY CONTROL METHOD

Process 2 levels of a quality control material (**Standatrol S-E 2 niveles**) with known uric acid concentrations, with each determination.

REFERENCE VALUES

Serum or plasma

Age	mg/dL	µmol/L
0-14 days	2.8 - 12.7	167 - 755
15 days - < 1 year	1.6 - 6.3	95 - 375
1 - 3 years	1.8 - 4.9	107 - 291
3 - 5 years	2.0 - 4.9	117 - 291
5 - 8 years	1.9 - 5.0	116 - 295
9 - 10 years	2.4 - 5.5	142 - 326
11-12 years	2.6 - 5.8	156 - 345
13-79 years (M)	3.7 - 7.7	218 - 459
13-79 years (F)	2.5 - 6.2	147 - 366

Urine

Average diet: 250750 mg/24 hours

Low purine diet:

Women: < 400 mg/24 hours

Men: < 480 mg/24 hours

High purine diet:

< 1000 mg/24 hours

It is recommended that each laboratory establish its own ranges or reference values, considering age, sex, dietary habits, and other factors.

UNIT CONVERSION

Uric acid (mg/dL) x 0.059 = Uric acid (mmol/L)

Uric acid (mg/24 hrs) x 0.0059 = Uric acid (mmol/24 hrs)

PROCEDURE LIMITATIONS

See known interfering substances in SAMPLE.

PERFORMANCE

The tests were performed on a CMD Series analyzer.

a) Precision:

It was tested using samples of human origin and following the protocol established by the CLSI EP5-A guide: for 20 days, 2 runs per day were performed. Each run was performed in duplicate (20x2x2; n = 80).

Sample	Mean (mg/dL)	Repeatability (CV _R %)	Intermediate precision (CV _{IL} %)
ORI	2.03	1.0%	1.2%
POOL	6.77	0.6%	1.0%
SE1	4.36	0.8%	1.3%
SE2	8.95	0.6%	1.0%

b) Detection capacity:

The Limits of Blank (LoB) and Detection (LoD) were tested according to a protocol based on EP17-A2 ($\alpha=\beta=0.05$). A daily quintuplicate run was performed for 5 days (1x5x5) at 6 sample levels with low uric acid concentration.

Parameter	Value (mg/dL)
LoB	0.01
LoD	0.02

c) Linearity:

Linearity testing was performed following the CLSI EP6-A protocol guidelines, following the weighted least squares (WLS) linear regression model with intercept. According to this model, the reaction is linear up to 20 mg/dL.

Measurement range is 0.02 mg/dL to 20 mg/dL.

d) Comparison of methods:

It was tested following the guidelines established by the CLSI EP9-A protocol. Assays with Uric Acid fast reagent were performed on the CMD800iX2, which was calibrated with **Calibrador A Plus** and controlled with **Standatrol S-E 2 niveles**. A similar method for uric acid determination was used as a comparative method. A total of 120 serum or plasma samples were run in duplicate. Passing-Bablok nonparametric regression was used for regression analysis.

$$y=0.988x+0.004$$

$$r=0.999$$

The sample concentration range was 2.01 to 19.40 mg/dL

PARAMETERS FOR AUTOMATED ANALYZERS

For programming instructions refer to the user's manual of the analyzer in use.

Request the applications for the analyzers marketed by Wiener lab. Applications not provided by Wiener lab must be validated.

WIENER LAB PROVIDES

- 2 x 50 mL w/ Standard (Code 1840112)
- 2 x 100 mL w/ Standard (Code 1840111)
- 4 x 100 mL w/ Standard (Code 1840113)
- 3 x 50 mL (Code 1008164)
- 4 x 50 mL (Code 1009403)
- 3 x 60 mL (Code 1009713)
- 4 x 60 mL (Code 1009298)

REFERENCES

- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, Sixth Edition. Elsevier 2018.
- Tietz Clinical Guide to Laboratory Tests, 4th edition. St. Louis (MO): Saunders Elsevier 2006.
- CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

- CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. Evaluation of Linearity of Quantitative Measurement Procedures. 2nd ed. CLSI guideline EP06. Clinical and Laboratory Standards Institute; 2020.
- CLSI. Interference Testing in Clinical Chemistry. 3rd ed. - CLSI guideline EP07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
- CLSI. Measurement Procedure Comparisons and Bias Estimation Using Patient Samples. 3rd ed. CLSI guideline EP09c. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.

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Wiener Laboratorios S.A.I.C.