Uricostat

enzimático

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Uricostat

enzimático

For uric acid determination in serum or plasma


SUMMARY

Uric acid is a metabolite found in purines, nucleic acids and nucleoproteins.

Serum uric acid concentration usually varies from one individual to another depending on several factors such as: sex, diet pattern, ethnic origin, genetic constitution, pregnancy. Abnormal levels of serum uric acid indicate metabolic disorders of the originating substances or inadequate excretion.

PRINCIPLE

The reaction scheme is as follows:

PROVIDED REAGENTS

A. Reagent A: vial containing more than 100 U peroxidase (POD), 12.5 mmol phosphates buffer for pH 7.3 and 3 mmol 4-aminophenazone.

B. Reagent B: 24 mmol/l chlorophenol solution.

C. Reagent C: uricase (UOD) solution greater than or equal to 3 kU/l.

S. Standard: 10 mg/dl uric acid solution.

Final concentrations

UOD.....	≥ 26 U/l
POD.....	≥ 200 U/l
Phosphates.....	25 mM
4-AP.....	6 mM
Chlorophenol.....	2.4 mM
pH.....	7.3 ± 0.1

NON-PROVIDED REAGENTS

Distilled water. See PROCEDURE LIMITATIONS.

INSTRUCTIONS FOR USE

Standard: ready to use.

Reagent A: reconstitute adding 50 ml distilled water, homogenize immediately by inversion to allow dissolution. Label and date. Before use, mix by inversion, without stirring.

Reagent B: ready to use.

Reagent C: ready to use.

Working Reagent: according to the work volume, place in a test tube 80 parts distilled water, 10 parts Reagent A (reconstituted), 10 parts Reagent B and 0.8 parts Reagent C. Mix by inversion without stirring, avoiding foam formation. Label and date. Different quantities may be prepared following the stated

proportions. Additionally, it is important to respect the order of reagent addition and ensure perfect homogenization, so that Reagent B is unable to deteriorate the Working Reagent.

WARNINGS

Reagents are for "in vitro" diagnostic use.

Reagent B: is harmful and corrosive. H301 + H311 + H331: Toxic if swallowed, in contact with skin or if inhaled. H314: Causes severe skin burns and eye damage. P262: Do not get in eyes, on skin, or on clothing. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302 + P352: IF ON SKIN: Wash with plenty of soap and water. P280: Wear protective gloves/protective clothing/eye protection/face protection.

Use the reagents according to the working procedures for clinical laboratories. The reagents and samples should be discarded according to the local regulations in force.

STABILITY AND STORAGE INSTRUCTIONS

Provided Reagents: stable in refrigerator (2-10°C) until the expiration date shown on the box. Do not expose to high temperatures for long periods.

Reagent A (reconstituted): stable for 3 months at room temperature (below 20°C) or for 10 months in refrigerator (2-10°C) from reconstitution date. At low temperatures this reagent may crystallize. In this case, before use, redissolve completely at room temperature and homogenize. The dissolution can be accelerated by placing the vial in water bath at 37°C for several minutes.

Working reagent: in a glass caramel colored bottle and in refrigerator (2-10°C) it is stable for 30 days from the date of preparation. At room temperature it is stable for 5 days.

INSTABILITY OR DETERIORATION OF REAGENTS

During use, the Working reagent may develop a slightly pink color, which does not affect the performance, provided that a Blank is run with each test batch and a Standard is run periodically.

Discard every time Blank readings are above 0.160 O.D. or when Standard readings are abnormally low.

SAMPLE

Serum or heparinized plasma

a) Collection: obtain serum or plasma as usual. Remove serum from clot as soon as possible within two hours from collection. If urine is used, it should be preferably fresh.

b) Known interfering substances:

- Icteric sera or with visible or intense hemolysis produce erroneous results, thus they should not be used.
- Fluoride inhibits uricase.
- Drugs: strong reducing substances such as ascorbic acid (vitamin C), buscapina (butyl bromide hyoscine), etc., supplied in large doses interfere in the reaction consuming H_2O_2 . For this reason, whenever possible, discontinue the medication to the patient 24 hours before sampling. Otherwise, falsely decreased values may be obtained.
- No interferences are observed by: bilirubin up to 120 mg/l, hyperlipidemia or slight hemolysis.

Refer to the bibliography of Young, D.S. in References for effect of drugs on the present method

d) Stability and storage instructions: use fresh serum. If the assay is not immediately performed, the serum may be stored up to 3 days in the freezer (without addition of preservatives).

REQUIRED MATERIAL (non-provided)

- Spectrophotometer or photocolorimeter.
- Photocolorimeter tubes or spectrophotometric square cuvettes.
- Caramel colored glass bottle.
- Water bath at 37°C (optional)
- Watch or timer.

ASSAY CONDITIONS

- Wavelength: 505 nm in spectrophotometer or in photocolorimeter with green filter (490-530 nm)
- Reaction temperature: 37°C or room temperature (18-25°C)
- Reaction time: 5 minutes at 37°C or 30 minutes at room temperature (18-25°C)
- Sample volume: 50 ul
- Working Reagent volume: 2.5 ml
- Final reaction volume: 2.55 ml

Sample and Reagent volumes may be proportionally decreased or increased (e.g. 20 ul Sample + 1 ml Working Reagent or 100 ul + 5 ml)

PROCEDURE

In three photocolorimetric tubes or spectrophotometer cuvettes labeled B (Blank), S (Standard) and U (Unknown), place:

	B	S	U
Standard	-	50 ul	-
Sample	-	-	50 ul
Working Reagent	2.5 ml	2.5 ml	2.5 ml

Mix gently and incubate in water bath at 37°C for 15 minutes or at room temperature (18-25°C) for 30 minutes. Remove, let cool and read in spectrophotometer with green filter (490-530 nm), setting the instrument to zero with the Blank.

CALCULATIONS

$$\text{Uric acid (mg/dl)} = U \times f \quad \text{where } f = \frac{10 \text{ mg/dl}}{S}$$

QUALITY CONTROL METHOD

Each time the test is performed, analyze two levels of a quality control material (**Standatrol S-E 2 niveles**) with known uric acid concentration.

REFERENCE VALUES

The following range is reported in normal adults with standard protein intake:

Men: 2.5-6.0 mg/dl

Women: 2.0-5.0 mg/dl

Uric acid levels vary according to age, sex and the different diets, and there has been observed a seasonal variation. It is therefore recommended that each laboratory establish its own reference intervals or values taking into account these factors.

UNITS CONVERSION

$$\text{Uric acid (mg/dl)} \times 0.059 = \text{Uric acid (mmol/l)}$$

PROCEDURE LIMITATIONS

See Known interfering substances under SAMPLE.

Reducing agents decrease the color response, while oxidizing agents color the Reagents increasing the Blanks. These agents are often found in the distilled water used to prepare the working reagent; therefore it is advisable to monitor its quality.

Detergents, heavy metals and cyanides are enzymatic inhibitors.

PERFORMANCE

a) Reproducibility: the following results were obtained processing replicates of the same samples in the same day:

Level	S.D.	C.V.
4.56 mg/dl	± 0.103 mg/dl	2.26 %
9.70 mg/dl	± 0.131 mg/dl	1.35 %

b) Recovery: by adding known amounts of uric acid to different sera, a recovery between 97 and 101% was obtained for 10 mg/dl uric acid level.

c) Detection limit: depending on the photometer used. In spectrophotometers, the 10 mg/dl Standard provides a reading of approximately 0.190 O.D., which means that the minimum detectable concentration change under these conditions for 0.001 O.D. will be approximately 0.05 mg/dl.

d) Linearity: the reaction is linear up to 20 mg/dl. For higher values, repeat the determination using half sample volume and multiplying the result by 2.

WIENER LAB. PROVIDES

Kit for 500 ml (Cat. N° 1840101).

REFERENCES

- Chu, S.Y. - Can. J. Med. Technol. 40/5:154 (1978).
- Day, J.H. - La Prensa Médica Argentina 58/15:786 (1971).

STABILITY OF FINAL REACTION

Final reaction color is stable for 45 minutes, thus absorbance should be read within that period.

- Donadon, V.; Barbieri, E.; Canterin, A. - LAB III/ 4:473 (1976).
- I.F.C.C. - Clin. Chim. Acta 87/3:459 F (1978).
- Lauriat, F. - Pharm. Biol. 11/108:163 (1977).
- Trinder, P. - Ann. Clin. Biochem. 6/24 (1969).
- Young, D.S. - "Effects of Drugs on Clinical Laboratory Tests", AACC Press, 4th ed., 2001.

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Wiener Laboratorios S.A.I.C.